Explanation of Symbols Which Appear on the Device Labeling
Explicación de símbolos que aparecen en la etiqueta del dispositivo
Explication des symboles qui apparaissent sur l'étiquette de l'appareil
Erklärung der Symbole auf dem Etikett dieses Gerätes
Spiegazione dei simboli che figurano nelle etichette del dispositivo

Caution or Warning
Precaución o advertencia
Avertissement ou mise en garde
Achtung oder Warnung
Attenzione o avvertenza

Refer to Instructions for Use
Consulte las instrucciones de uso.
Se reporter au mode d'emploi
Siehe Gebrauchsanleitung
Consultare le Istruzioni per l'uso

Manufacturer
Fabricante
Fabricant
Hersteller
Fabbricante

CE Certification Mark
Marca de certificación CE
Marquage de certification CE
CE Zertifizierungsmarke
Marcatura CE

Part number
Número de pieza
Référéncie
Teilenummer
Numero di articolo

Lot number
Número de lote
Numéro du lot
Losnummer
Numero di lotto

This Product is not manufactured with Natural Rubber Latex.
Este producto no se fabrica con látex de caucho.
Ce produit est fabriqué sans latex
Dieses Produkt ist nicht mit Kautschuklatex hergestellt
Il prodotto non contiene lattice.

Date of Manufacture
Fecha de fabricación
Date de fabrication
Herstellungsdatum
Data di produzione

EU Rep: Medical Product Services GmbH
Borngasse 20, 35619 Braunfels Germany
info@mps-gmbh.eu

Innova® Full Face Vented Mask Instruction Guide

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






**These instructions contain the following languages:
English, Spanish, French, German, and Italian**

Sleepnet Corporation
5 Merrill Industrial Drive
Hampton, NH 03842 USA
Tel: (603) 758-6600
Fax: (603) 758-6699
Toll Free: 1-800-742-3646 (USA AND CANADA ONLY)
www.sleepnetmasks.com

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

en

INNOVA FULL FACE VENTED MASK INSTRUCTION GUIDE
 CONTENTS: INNOVA FULL FACE VENTED MASK WITH HEADGEAR
 MADE IN THE U.S.A. OF U.S. AND IMPORTED PARTS

RX ONLY
INTENDED USE

The Innova Full Face Vented Mask is intended to be used with positive airway pressure devices, such as CPAP or bilevel operating at or above 3 cm H₂O. The masks are to be used on adult patients (>30kg) for whom positive airway pressure therapy has been prescribed. The mask is intended for single patient use in the home, hospital or institutional environment.

NOTES

- The mask system does not contain PVC, DEHP, or phthalate materials. If the patient has ANY reaction to any part of the mask system, discontinue use.
- This product is not manufactured with natural rubber latex.
- An exhalation port is built into this mask so a separate exhalation port is not required.
- The mask includes an anti-asphyxia valve to help reduce rebreathing of exhaled air in the event that the positive pressure device does not work properly.

CAUTIONS

- In the U.S. Federal Law restricts this device for sale by or on the order of a physician.
- Do not expose the mask to temperatures above 140°F (60°C).
- The mask is intended for SINGLE PATIENT USE ONLY. It may be cleaned and used repeatedly on the same person, but must not be used on multiple persons. Do not sterilize or disinfect.

WARNINGS

- This mask is not suitable for providing life support ventilation.
- At low pressures, the flow through the exhalation ports may be inadequate to clear all exhaled gas from the tubing. Some rebreathing may occur.
- Do not use this mask with baseline pressure less than 3 cm H₂O.
- This mask should be used with the positive airway pressure device recommended by a physician or respiratory therapist.
- Do not block or try to seal the vent holes or non-rebreathing valve. The vent holes allow a continuous flow of air out of the mask.
- The mask should not be used unless the positive airway pressure device is turned on and operating properly. Explanation of the warning: When the positive airway pressure device is turned on and functioning properly, new air from the device flushes the exhaled air out through the mask's vent holes. When the positive airway pressure device is not operating, exhaled air may be rebreathed. Rebreathing of exhaled air can, in some circumstances, lead to suffocation. This warning applies to most models of CPAP machines and masks.
- If oxygen is used with this CPAP mask, the oxygen flow must be turned off when the CPAP machine is not operating. Explanation of the warning: When the CPAP device is not in operation, and the oxygen flow is left on, oxygen delivered into the device tubing may accumulate within the CPAP machine enclosure. Oxygen accumulated in the CPAP machine enclosure will create a risk of fire.

Note: At a fixed flow rate of supplemental oxygen, the inhaled oxygen concentration will vary, depending on pressure settings, patient breathing, and leak rate.

- No smoking or open flames such as candles may be used when oxygen is in use.
- To minimize the risk of vomiting during sleep, avoid eating or drinking for three (3) hours before using the mask.
- Any unusual skin irritation, chest discomfort, shortness of breath, gastric distension, abdominal pain, belching or flatulence from ingested air, or severe headache experienced during or immediately after use should be reported to a healthcare professional.
- Using a mask may cause tooth, gum or jaw soreness or aggravate an existing dental condition. If symptoms occur, consult a physician or dentist.
- Consult a physician if you experience any of the following symptoms while using the mask: Drying of the eyes, eye pain, eye infections, or blurred vision. Consult an ophthalmologist if symptoms persist.
- Do not overtighten the headgear straps. Signs of overtightening include: excessive redness, sores, or bulging skin around the edges of the mask. Loosen headgear straps to alleviate symptoms.
- Significant unintentional leak may occur without properly fitting the mask.

CONTRAINDICATIONS

Do not use this product if you are experiencing nausea, vomiting, taking a prescription drug that may cause vomiting, or if you are unable to remove the mask on your own.

BEFORE USE

- Inspect mask and valve prior to each use. Replace the mask if parts are damaged or gel is exposed

due to tear or punctures.

- ▶ Check the anti-asphyxia valve. With the positive pressure device turned off, verify that the valve flapper is positioned so that room air can flow through the large openings in the valve (Figure 1a). Turn on the system; the flapper should close, and system air should flow through the mask (Figure 1b). If the flapper does not close or does not function properly, replace the mask. Do not block the opening on the anti-asphyxia valve or the vent holes. Ensure that the valve is not blocked with secretions and that the flapper is dry.

GENERAL INFORMATION

This mask has a non-transferable three month limited warranty for manufacturer's defects from the date of purchase by the initial consumer. If a mask used under normal conditions should fail, Sleepnet will replace the mask or its components. For additional information about the Sleepnet Warranty, please visit <http://www.sleepnetmasks.com/products/our-promise.cfm>.

NOTE: Please dispose of the mask and accessories in accordance with local and regional guidelines.

CLEANING AND MAINTENANCE-SINGLE PATIENT USE

Cleaning is recommended on a daily basis. Hand wash the mask and headgear (as needed) in warm water using a mild detergent. Rinse thoroughly after washing. Visually verify that the vent holes are clear and allow the mask to air dry. Washing your face prior to fitting the mask will remove excess facial oils and help extend the life of the mask.

Weekly Cleaning (if desired): The Innova Mask is dishwasher safe if placed in top rack; small components should be secured. Air dry before use.

NOTE: The gel cushion, although silky to the touch, will not tear, leak, or detach from the mask under normal use and following recommended cleaning and maintenance protocols. However, the gel cushion is not indestructible. It will tear, cut, and rip if subjected to mishandling or abuse. The warranty does not cover these types of product failures. Therefore, please handle your new mask with care.

SPECIFICATIONS

Resistance to air flow at 50 L/min: 0.07

Resistance to air flow at 100 L/min: 0.39

DEAD SPACE VOLUME (mL approx.): Small - 155.3 ml/ Medium - 181.7ml/ Large - 196.6ml

Open-to-atmosphere pressure: 2.1 cmH₂O

Closed-to-atmosphere pressure: 2.28 cmH₂O

Operating Temperature: Do not expose the mask to temperatures above 140°F (60°C).

EXPECTED SERVICE LIFE

- ▶ Mask Frame Components - 6 months
- ▶ Headgear - 6 months

USER INSTRUCTIONS

1. Remove the mask and the headgear from the package.
2. Place the headgear on a table, label side down with the narrowest part close to you.
3. Set the mask with the gel side down onto the center of the headgear.
4. Take each headgear strap and connect it to its respective location on the mask by looping the Velcro strap through the headgear slots and attaching it back to itself.
5. Once the headgear has been properly attached to the mask, follow the fitting tips below to fit the mask:
6. Detach one Velcro® strap by sliding the Velcro loop off the quick release hook, slip the headgear over the head and center the mask on the face (Figure 2).
7. Find the detached Velcro® strap and slide the Velcro loop back on the quick release hook (Figure 3).
8. Be sure the headgear straps are parallel and the headgear is fully extended and flat around the crown of the head (Figure 4).
9. Adjust the bottom headgear straps by using the Velcro® tabs. Make sure the Active Headgear Connector is parallel to the face and that there is equal tension from the headgear on both sides. After tightening, make sure the Active Headgear Connector™ can move freely. (Figure 5).

Note: Do NOT overtighten the headgear.

10. Adjust the upper headgear straps by using the Velcro® tabs (Figure 6).

11. Adjust the bottom headgear straps to make sure the bottom of the mask is seated comfortably on the face.

Note: For optimal seal and comfort, mask should be snug, but not tight. Adjust the tension by detaching the Velcro straps, tightening or loosening the straps, and reattaching the Velcro.

12. Connect the mask tubing assembly to the positive airway pressure device tubing and turn the device on to the prescribed pressure (Figure 7).
13. Turn head from side to side while lying on the bed and identify areas of leak that may prevent sleep.

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Velcro® is a registered trademark of Velcro Industries B.V.
Innova® is a trademark of Sleepnet Corporation.
CE₀₁₂₃ Velcro® est une marque déposée de Velcro Industries B.V.
Innova® est une marque commerciale de Sleepnet Corporation.