Explanation of Symbols Which Appear on the Device Labeling  
Explicación de los símbolos que aparecen en la etiqueta del dispositivo  
Explication des symboles qui apparaissent sur l'étiquetage de l'appareil  
Пояснение символов, используемых для маркировки устройства

Caution or Warning  
Precaución o advertencia  
Mise en garde ou Avertissement  
Внимание или Осторожно

Refer to Instructions for Use  
Consulte las instrucciones de uso  
Consulter les instructions d'utilisation  
См. инструкции по применению

CE Certification Mark  
Certificación CE  
Marquage de certification CE  
Знак сертификации CE

Lot number  
Número del lote  
Numéro de lot  
Номер партии

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Contents: Ascend Full Face Vented Mask, replaceable cushions, and headgear

Made in the U.S.A. of U.S. and imported parts

Rx only

### INTENDED USE

The Ascend Full Face Vented Mask is intended to be used with positive airway pressure devices, such as CPAP or bilevel operating at or above 3 cm H<sub>2</sub>O. The mask is to be used on adult patients (>30kg) for whom positive airway pressure therapy has been prescribed. The mask is intended for single patient use in the home, hospital or institutional environment.

### NOTES

- The mask system does not contain PVC, DEHP, or phthalate materials. If the patient has ANY reaction to any part of the mask system, discontinue use.
- This product is not manufactured with natural rubber latex.
- An exhalation port is built into this mask so a separate exhalation port is not required.
- The mask includes an anti-asphyxia valve to help reduce rebreathing of exhaled air in the event that the positive pressure device does not work properly.

### CAUTIONS

- In the U.S. Federal Law restricts this device for sale by or on the order of a physician.
- Do not expose the mask to temperatures above 140°F (60°C).
- The mask is intended for SINGLE PATIENT USE ONLY. It may be cleaned and used repeatedly on the same person, but must not be used on multiple persons. Do not sterilize or disinfect.

### WARNINGS

- This mask is not suitable for providing life support ventilation.
- At low pressures, the flow through the exhalation ports may be inadequate to clear all exhaled gas from the tubing. Some rebreathing may occur.
- Do not use this mask with baseline pressure less than 3 cm H<sub>2</sub>O.
- This mask should be used with the positive airway pressure device recommended by a physician or respiratory therapist.
- Do not block or try to seal the vent holes or non-rebreathing valve. The vent holes allow a continuous flow of air out of the mask.
- The mask should not be used unless the positive airway pressure device is turned on and operating properly. Explanation of the warning: When the positive airway pressure device is turned on and functioning properly, new air from the device flushes the exhaled air out through the mask's vent holes. When the positive airway pressure device is not operating, exhaled air may be rebreathed. Rebreathing of exhaled air can, in some circumstances, lead to suffocation. This warning applies to most models of CPAP machines and masks.
- If oxygen is used with this CPAP mask, the oxygen flow must be turned off when the CPAP machine is not operating. Explanation of the warning: When the CPAP device is not in operation, and the oxygen flow is left on, oxygen delivered into the device tubing may accumulate within the CPAP machine enclosure. Oxygen accumulated in the CPAP machine enclosure will create a risk of fire.

Note: At a fixed flow rate of supplemental oxygen, the inhaled oxygen concentration will vary, depending on pressure settings, patient breathing, and leak rate.

- No smoking or open flames such as candles may be used when oxygen is in use.
- To minimize the risk of vomiting during sleep, avoid eating or drinking for three (3) hours before using the mask.
- Any unusual skin irritation, chest discomfort, shortness of breath, gastric distension, abdominal pain, belching or flatulence from ingested air, or severe headache experienced during or immediately after use should be reported to a healthcare professional.
- Using a mask may cause tooth, gum or jaw soreness or aggravate an existing dental condition. If symptoms occur, consult a physician or dentist.
- Consult a physician if you experience any of the following symptoms while using the mask: Drying of the eyes, eye pain, eye infections, or blurred vision. Consult an ophthalmologist if symptoms persist.

- Do not overtighten the headgear straps. Signs of overtightening include: excessive redness, sores, or bulging skin around the edges of the mask. Loosen headgear straps to alleviate symptoms.
- Significant unintentional leak may occur without properly fitting the mask or attaching the mask cushion.

#### CONTRAINDICATIONS

Do not use this product if you are experiencing nausea, vomiting, taking a prescription drug that may cause vomiting, or if you are unable to remove the mask on your own.

#### BEFORE USE

- Inspect mask and valve prior to each use. Replace the mask if parts are damaged or gel is exposed due to tear or punctures.
- Check the anti-asphyxia valve. With the positive pressure device turned off, verify that the valve flapper is positioned so that room air can flow through the large openings in the valve (Figure 1a). Turn on the system; the flapper should close, and system air should flow through the mask (Figure 1b). If the flapper does not close or does not function properly, replace the mask. Do not block the opening on the anti-asphyxia valve or the vent holes. Ensure that the valve is not blocked with secretions and that the flapper is dry.

#### GENERAL INFORMATION

This mask shell has a non-transferable three month limited warranty (or maximum service life, if shorter) for manufacturer's defects from the date of purchase by the initial consumer. If a mask used under normal conditions should fail, Sleepnet will replace the mask or its components. For additional information about the Sleepnet Warranty, please visit [www.sleepnetmasks.com](http://www.sleepnetmasks.com). Per EU MDR 2017/745, please report serious incidents to manufacturer and EU Member State competent authority.

NOTE: The mask shell and headgear have a service life of 6 months. The cushion has a service life of 1 month.

NOTE: Please dispose of the mask and accessories in accordance with local and regional guidelines.

#### TECHNICAL INFORMATION

Resistance to air flow at 50 L/min: 0.58 cm H<sub>2</sub>O

Resistance to air flow at 100 L/min: 1.02 cm H<sub>2</sub>O

DEAD SPACE VOLUME (mL approx.): Small 165 Medium 212 Large 221

Open-to-atmosphere pressure: 1.92 cm H<sub>2</sub>O

Closed-to-atmosphere pressure: 1.85 cm H<sub>2</sub>O

Operating Temperature: Do not expose the mask to temperatures above 140°F (60°C).

#### ASSEMBLY

1. Remove the mask and the headgear from the package.
2. The Ascend Full Face Mask is packaged with three cushion sizes. The medium size has been assembled. Choose the smallest size that is wide enough to extend beyond the nostrils and lips but does not obstruct breathing (Figure 2).
3. If a new size is needed, separate the cushion from the mask frame by gently squeezing the cushion ring (not the AIRgel cushion) near the left and right alignment pins while pulling away from the mask (Figure 3).
4. Ensure the new cushion is oriented correctly by aligning all three alignment pins with the notches located in the three corners of the mask (Figure 4)
5. Gently push the cushion alignment pins into the notches. Make sure that the Support Tabs stay outside the Mask Shell.
6. Inspect to make sure that each pin is pushed completely into the notch (Figure 4).
7. Place the headgear on a table, label side down with the narrowest part close to you.
8. Set the mask with the gel side down onto the center of the headgear.
9. Take each headgear strap and connect it to its respective location on the mask by looping the Velcro® strap through the headgear slots and attaching it back to itself (Figure 5).
10. Once the headgear and cushion have been properly attached to the mask, follow the fitting tips below to fit the mask.

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## FITTING

- Detach one Velcro® strap by sliding the Velcro loop off the quick release hook (Figure 6), seat the mask at the center of the chin and tip the mask up to contact the nose. The mask may sit lower on the nose than expected, depending on facial features (Figure 7). Slip the headgear over the head and center the mask on the face.
- Find the detached Velcro® strap and slide the Velcro loop back on the quick release hook (Figure 8).
- Adjust the bottom headgear straps by using the tabs. After tightening, make sure the Active Headgear Connector can move freely (Figure 9).  
NOTE: Do NOT overtighten the headgear.
- Adjust the upper headgear straps by using the tabs (Figure 10).  
NOTE: The Flexible Spacebar should not touch the forehead (Figure 11).
- Be sure the headgear straps are parallel and the headgear is fully extended and flat across the back of the head. (Figure 8)  
NOTE: For optimal seal and comfort, the mask should be snug, but not tight. Adjust the tension by detaching the Velcro® straps, tightening or loosening the straps, and reattaching the Velcro.
- Connect the mask to the positive airway pressure device tubing and breathe normally, while moving head from side-to-side (Figure 12).
- If there are leaks, adjust the headgear straps closest to the area of leak.
- If leaks persist, try moving the mask higher or lower around the nose to find the best fit.

## CLEANING AND MAINTENANCE

Cleaning is recommended on a daily basis. Disassemble the mask cushion according to the disassembly instructions. Hand wash the mask and components in warm water using a mild detergent. Rinse the mask thoroughly after washing. Visually verify that the vent holes are clear and allow the mask to air dry. When the mask and cushion are dry, reassemble according to the reassembly instructions. Washing your face prior to fitting the mask will remove excess facial oils and help extend the life of the mask.

HEADGEAR: Hand wash the headgear in warm water using mild detergent. Rinse after washing and allow to air dry.

NOTE: Inspect mask daily or prior to each use. Replace the mask if vent holes are blocked, if parts are damaged, or if gel cushion has tears or punctures.

NOTE: The gel cushion, although silky to the touch, will not tear, leak, or detach from the mask under normal use. The gel cushion is not indestructible. It will tear, cut, and rip if subjected to mishandling or abuse.

Velcro® is a registered trademark of Velcro Industries B.V.

Ascend® is a registered trademark of Sleepnet Corporation.

Visit our website at [www.sleepnetmasks.com](http://www.sleepnetmasks.com).

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