Explanation of Symbols Which Appear on the Device Labeling
 Explicación de los símbolos que aparecen en la etiqueta del dispositivo
 Erklärung der Symbole auf dem Etikett dieses Gerätes
 Explication des symboles qui apparaissent sur l'étiquette de l'appareil
 Spiegazione dei simboli che figurano sulle etichette del dispositivo
 Laitteen etiketissä olevien symbolien selitykset
 装置标签上的符号说明
 装置のラベルに表示される記号の説明
 Cihaz Etiketinde Görünen Sembollerin Açıklaması
 Пояснение символов, используемых для маркировки устройства
 Explicații cu privire la simbolurile ce apar pe eticheta dispozitivului

Veraseal® 2 Full Face Vented Mask

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Australian Sponsor:
 Breathing and Medical.
 127 Bowden St, Meadowbank, NSW 2114
 www.breathingandmedical.com.au
 1800 335 333

Medical Product Services GmbH
 Borngasse 20, 35619 Braunfels Germany
 info@mps-gmbh.eu

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
These instructions contain the following languages:

English, Spanish, German, French, Italian, Finnish, Chinese, Japanese, Turkish, Russian, and Romanian


Sleepnet Corporation
 5 Merrill Industrial Drive
 Hampton, NH 03842 USA
 Tel: (603) 758-6600
 Fax: (603) 758-6699
 Toll Free: 1-800-742-3646 (USA)
 www.sleepnetmasks.com


1a


1b

2


4


3


5

6



6a

This page intentionally left blank.

7

This page intentionally left blank.

Veraseal[®] 2 Vented Mask

1. Flexible Spacebar™
2. Mask Shell with vent holes
3. Airgel[®] Cushion
4. Active Headgear Connector
5. Elbow Assembly with AAV
6. Oxygen port
7. Quick Release Hook

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Contents: Veraseal 2 Vented Mask with headgear

Made in the U.S.A. of U.S. and imported parts

Rx only

INTENDED USE

The Veraseal® 2 Full Face Vented Mask is a disposable mask intended to be used with positive airway pressure devices, such as CPAP or bi-level, operating at or above 3 cm H₂O. The mask is to be used on adult patients (>30kg) for whom positive airway pressure therapy has been prescribed. The mask is intended for **SHORT TERM SINGLE PATIENT USE** (maximum of 7 days) in the hospital or institutional environment.

NOTES

- The mask system does not contain PVC, DEHP, or phthalates. If the patient has ANY reaction to any part of the mask system, discontinue use.
- This product is not manufactured with natural rubber latex.
- An exhalation port is built into this mask so a separate exhalation port is not required.
- The mask includes an anti-asphyxia valve to help reduce rebreathing of exhaled air in the event that the positive pressure device does not work properly.
- Verify that the mask is the right size.

CONTRAINDICATION

This product should not be used if the patient is experiencing nausea, vomiting, taking a prescription drug that may cause vomiting, or if they are unable to remove the mask on their own.

CAUTIONS

- In the U.S. Federal Law restricts this device for sale by or on the order of a physician.
- Do not expose the mask to temperatures above 140°F (60°C).
- The mask is intended for **SINGLE PATIENT USE ONLY**. It is intended for short term use (not more than 7 days) and must not be used on multiple persons. Do not sterilize or disinfect.
- Consult attending physician prior to using mask if the patient is using any medications or devices to remove deep secretions.

WARNINGS

- This mask is not suitable for providing life support ventilation.
- The elbow and anti-asphyxia valve assembly have specific safety functions. The mask should not be used if the valve is damaged. The mask should be discarded if the valve is damaged, distorted, or torn.
- If the patient has respiratory insufficiency, they should wear this mask **ONLY** when therapy is being delivered.
- At low pressures, the flow through the exhalation ports may be inadequate to clear all exhaled gas from the tubing. Some rebreathing may occur.
- Use this mask with baseline pressure of 3 cm H₂O or higher.
- Do not block or try to seal the vent holes or anti-asphyxia valve. The vent holes allow a continuous flow of air out of the mask.
- The mask should not be used unless the positive airway pressure device is turned on and operating properly. Explanation of the warning: When the positive airway pressure device is turned on and functioning properly, new air from the device flushes the exhaled air out through the mask's vent holes. When the positive airway pressure device is not operating, exhaled air may be rebreathed. Rebreathing of exhaled air can, in some circumstances, lead to suffocation. This warning applies to most models of CPAP machines and masks.
- If oxygen is used with this CPAP mask, the oxygen flow must be turned off when the CPAP machine is not operating. Explanation of the warning: When the CPAP device is not in operation, and the oxygen flow is left on, oxygen delivered into the device tubing may accumulate within the CPAP machine enclosure. Oxygen accumulated in the CPAP machine enclosure will create a risk of fire.
- At a fixed flow rate of supplemental oxygen, the inhaled oxygen concentration will vary, depending on pressure settings, patient breathing, and leak rate.
- No smoking or open flames such as candles may be used when oxygen is in use.
- This mask is not intended for patients without a spontaneous respiratory drive. This mask should not be used on patients who are uncooperative, obtunded, unresponsive, or unable to remove the mask. This mask may not be suitable for persons with the following conditions: impaired cardiac sphincter function, excessive reflux, impaired cough reflex, and hiatal hernia.
- Any unusual skin irritation, chest discomfort, shortness of breath, gastric distension, abdominal pain, belching or flatulence from ingested air, or severe headache experienced by the patient during or immediately after use should be immediately reported to the attending healthcare professional.
- Using a mask may cause tooth, gum or jaw soreness or aggravate an existing dental condition. If symptoms occur, consult a physician or dentist.
- Consult a physician if the patient experiences any of the following symptoms while using the mask: Drying of the eyes, eye pain, eye infections, or blurred vision. Consult an ophthalmologist if symptoms persist.
- Do not overtighten the headgear straps. Signs of overtightening include: excessive redness, sores, or bulging skin around the edges of the mask. Loosen headgear straps to alleviate symptoms.
- Significant unintentional leak may occur without properly fitting the mask.

GENERAL INFORMATION

The mask has a service life of 7 days. Please dispose of the mask and accessories in accordance with local and regional guidelines. This mask and all its accessories do not contain natural rubber latex. Per EU MDR 2017/745, please report serious incidents to manufacturer and EU Member State competent authority.

CLEANING AND MAINTENANCE

This mask is disposable and is intended for single patient short term use (maximum of 7 days) only. The mask should be discarded if heavily soiled. The mask may be cleaned daily or as needed using a 70% v/v isopropyl alcohol swab or by wiping with a damp, soapy cloth. If using a damp cloth, use mild soap. Rinse thoroughly and allow to air dry, out of direct sunlight.

NOTE: The gel cushion, although silky to the touch, will not tear, leak, or detach from the mask under normal use. The gel cushion is not indestructible. It will tear, cut, and rip if subjected to mishandling or abuse.

TECHNICAL INFORMATION

Resistance to air flow at 50 L/min: 0.07 cm H₂O

Resistance to air flow at 100 L/min: 0.39 cm H₂O

DEAD SPACE (approx.): **Small** 155.3 ml **Medium** 181.7 ml **Large** 196.6 ml **Extra Large** 200.6 ml

Anti-asphyxia valve: Open-to-atmosphere pressure: 2.1 cm H₂O

Closed-to-atmosphere pressure: 2.28 cm H₂O

Operating Temperature: Do not expose the mask to temperatures above 140°F (60°C).

BEFORE USE

- ▶ Inspect mask daily or prior to each use. Replace the mask if vent holes are blocked, if parts are damaged, or if gel is exposed due to tears or punctures.
- ▶ Check the anti-asphyxia valve. With the positive pressure device turned off, verify that the valve flapper is positioned so that room air can flow through the large opening in the valve (Figure 1a). Turn on the system; the flapper should close, and system air should flow through the mask (Figure 1b). If the flapper does not close or does not function properly, replace the mask. Do not block the opening on the anti-asphyxia valve or the vent holes. Ensure that the valve is not blocked with secretions and that the flapper is dry.

MASK FITTING INSTRUCTIONS

1. Detach one Velcro® strap by sliding the Velcro loop off the quick release hook, slip the headgear over the head and center the mask on the face (Figure 2).
2. Find the detached Velcro® strap and slide the Velcro loop back on the quick release hook (Figure 3).
3. Be sure the headgear straps are parallel and the headgear is fully extended and flat across the back of the head (Figure 4).
4. Adjust the bottom headgear straps by using the Velcro® tabs. After tightening, make sure the Active Headgear Connector™ can move freely, (Figure 5).
5. Adjust the upper headgear straps by using the Velcro® tabs (Figure 6).

NOTE: Do NOT overtighten the headgear. The Flexible Spacebar should not touch the forehead (Figure 6a).

6. Adjust the bottom headgear straps to make sure the bottom of the mask is sealed comfortably against the patient's face.

NOTE: For optimal seal and comfort, mask should be snug, but not tight. Adjust the tension by detaching the Velcro straps, tightening or loosening the straps, and reattaching the Velcro.

7. Connect the mask tubing assembly to the positive airway pressure device tubing and turn the device on to the prescribed pressure (Figure 7).

Velcro® is a registered trademark of Velcro Industries B.V.

Veraseal® is a trademark of Sleepnet Corporation.

Visit our website at www.sleepnetmasks.com.

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