Explanation of Symbols Which Appear on the Device Labeling
Explicación de los símbolos que aparecen en la etiqueta del dispositivo
Erklärung der Symbole auf der Gerätezeichnung
Explication des symboles qui apparaissent sur l'étiquette de l'appareil
Illustrazione dei simboli visualizzati nell'etichetta del dispositivo
Laitteen merkinnoissa näkyvien symbolien selitys
Cihaz Etiketinde Görünen Sembollerin Açıklaması
Explicații cu privire la simbolurile ce apar pe eticheta dispozitivului

Caution or Warning
Precaución o advertencia
Vorsicht und Warnhinweis
Avertissement ou mise en garde
Precauzioni o Avvertenze
Huomio tai varoitus
İkaz veya Uyarı
Precauții sau avertismente

Refer to Instructions for Use
Consulte las instrucciones de uso.
Siehe die Gebrauchsanleitung
Se reporter au Mode d'emploi
Fare riferimento alle istruzioni d'uso
Katsa käyttöohjeita
Kullanım Talimatlarına Bakın
Consultați instrucțiunile de utilizare

CE Certification Mark
Marca de certificación CE
CE-Kennzeichnung
Marquage CE
Marchio di Certificazione CE
CE-sertifiointimerkintä
CE Belgelandirne İşareti
Marcă de certificare CE

Lot number
Número de lote
Chargennummer
Numéro du lot
Numero del lotto
Eränumero
Parti numarası
Număr lot

Manufacturer
Fabricante
Hersteller
Fabricant
Fabbricante
Valmistaja
Üretici
Producător

Part number
Referencia
Teilenummer
Référence
Numero di articolo
Osanumero
Parça numarası
Număr piesă

Medical Device
Producto sanitario
Medizinisches Gerät
Dispositif médical
Dispositivo medico
Lääkinnällinen laite
Medikal Cihaz
Dispozitiv medical

Date of Manufacture
Fecha de fabricación
Herstellungsdatum
Date de fabrication
Data di fabbricazione
Valmistuspäivämäärä
Üretim Tarihi
Data fabricației

Use by Date
Fecha de caducidad
Verfallsdatum
Date limite d'utilisation
Da utilizzare entro la data indicata
Viimeinen käyttöpäivä
Son Kullanma Tarihi
A se folosi până la data

This Product is not manufactured with Natural Rubber Latex.
Este producto no ha sido fabricado con látex de caucho natural.
Dieses Produkt wurde nicht mit Naturkautschukgummi hergestellt.
Ce produit est fabriqué sans latex de caoutchouc naturel.
Questo prodotto non è fabbricato in lattice di gomma naturale.
Tätä tuotetta ei ole valmistettu luonnonkumilateksista.
Bu Ürün, Doğal Kauçuk Lateks ile üretilmemiştir
Acest produs nu este fabricat cu latex din cauciuc natural.

Medical Product Services GmbH
Borngasse 20, 35619 Braunsfels
Germany
info@mps-gmbh.eu

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These instructions contain the following languages:

English, Spanish, German, French, Italian, Finnish, Turkish, Romanian

Sleepnet Corporation
5 Merrill Industrial Drive
Hampton, NH 03842 USA
Tel: (603) 758-6600
Fax: (603) 758-6699
Toll Free: 1-800-742-3646 (USA ONLY)
www.sleepnetmasks.com


Figure 1


Figure 2


Figure 3


Figure 4

This page intentionally left blank.


Figure 5


Figure 6


Figure 7a



Figure 7b

This page intentionally left blank.

Veraseal® 3 Vented Mask

1. Forehead Pad
2. Adjustment Screw
3. Vent Holes
4. Airgel® Cushion
5. Elbow Assembly with AAV
6. Oxygen port
7. Headgear Connector
8. Quick Release Hook
9. Headgear

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INTENDED USE

The Veraseal® 3 Full Face Vented Mask is a disposable mask intended to be used with positive airway pressure devices, such as CPAP or bi-level, operating at or above 3 cm H₂O. The mask is to be used on adult patients (>30kg) for whom positive airway pressure therapy has been prescribed. The mask is intended for SHORT TERM SINGLE PATIENT USE (maximum of 7 days) in the hospital or institutional environment.

NOTES

- The mask system does not contain PVC, DEHP, or phthalates.
- This product is not manufactured with natural rubber latex.
- If the patient has ANY reaction to any part of the mask system, discontinue use.
- An exhalation port is built into this mask so a separate exhalation port is not required.
- The mask includes an anti-asphyxia valve (air entrainment valve) to help reduce rebreathing of exhaled air in the event that the positive pressure device does not work properly.
- Verify that the mask is the right size by using the sizing guide.

CONTRAINDICATIONS

This product should not be used if the patient is experiencing nausea, vomiting, taking a prescription drug that may cause vomiting, or if they are unable to remove the mask on their own.

CAUTIONS

- In the U.S. Federal Law restricts this device for sale by or on the order of a physician.
- Do not expose the mask to temperatures above 140°F (60°C).
- The mask is intended for SINGLE PATIENT USE ONLY. It is intended for short term use (not more than 7 days) and must not be used on multiple persons. Do not sterilize or disinfect.
- Consult attending physician prior to using mask if the patient is using any medications or devices to remove deep secretions.

WARNINGS

- This mask is not suitable for providing life support ventilation.
- If the patient has respiratory insufficiency, they should wear this mask ONLY when therapy is being delivered.
- Use this mask with baseline pressure of 3 cm H₂O or higher.
- If oxygen is used with this CPAP mask, the oxygen flow must be turned off when the CPAP machine is not operating. Explanation of the warning: When the CPAP device is not in operation, and the oxygen flow is left on, oxygen delivered into the device tubing may accumulate within the CPAP machine enclosure. Oxygen accumulated in the CPAP machine enclosure will create a risk of fire.
- At a fixed flow rate of supplemental oxygen, the inhaled oxygen concentration will vary, depending on pressure settings, patient breathing, and leak rate.
- No smoking or open flames may be used when oxygen is in use.
- This mask is not intended for patients without a spontaneous respiratory drive. This mask should not be used on patients who are uncooperative, obtunded, unresponsive, or unable to remove the mask. This mask may not be suitable for persons with the following conditions: impaired cardiac sphincter function, excessive reflux, impaired cough reflex, and hiatal hernia.
- Any unusual skin irritation, chest discomfort, shortness of breath, gastric distension, abdominal pain, belching or flatulence from ingested air, or severe headache experienced by the patient during or immediately after use should be immediately reported to the attending healthcare professional.
- Using a mask may cause tooth, gum or jaw soreness or aggravate an existing dental condition. If symptoms occur, consult a physician or dentist.
- Consult a physician if the patient experiences any of the following symptoms occur while using the mask: drying of the eyes, eye pain, eye infections, or blurred vision. Consult an ophthalmologist if symptoms persist.
- Do not overtighten the headgear straps. Signs of overtightening include: excessive redness, sores, or bulging skin around the edges of the mask. Loosen headgear straps to alleviate symptoms.
- Significant unintentional leak may occur without properly fitting the mask. Follow the MASK FITTING INSTRUCTIONS provided.
- The elbow and anti-asphyxia valve assembly have specific safety functions. The mask should not be used if the valve flapper is damaged, distorted, or torn.

- At low pressures, the flow through the exhalation ports may be inadequate to clear all exhaled gas from the tubing. Some rebreathing may occur.
- Do not block or try to seal the exhalation port (vent holes) or anti-asphyxia valve (air entrainment valve). The vent holes allow a continuous flow of air out of the mask.
- The mask should not be used unless the positive airway pressure device is turned on and operating properly. When the positive airway pressure device is turned on and functioning properly, new air from the device flushes the exhaled air out through the mask's vent holes. When the positive airway pressure device is not operating, exhaled air may be rebreathed. Rebreathing of exhaled air can, in some circumstances, lead to suffocation. This warning applies to most models of CPAP machines and masks.

GENERAL INFORMATION

The mask has a service life of 7 days. Please dispose of the mask and accessories in accordance with local and regional guidelines. Per EU MDR 2017/745, please report serious incidents to manufacturer and EU Member State competent authority.

CLEANING AND MAINTENANCE

This mask is disposable and is intended for single patient short term use (maximum of 7 days) only. The mask should be discarded if heavily soiled. The mask may be cleaned daily or as needed using a 70% isopropyl alcohol swab or by wiping with a damp, soapy cloth. If using a damp cloth, use mild soap such as Ivory®. Rinse thoroughly and allow to air dry, out of direct sunlight.

NOTE: The gel cushion, although silky to the touch, will not tear, leak, or detach from the mask under normal use. The gel cushion is not indestructible. It will tear, cut, and rip if subjected to mishandling or abuse.

TECHNICAL INFORMATION

Resistance to air flow at 50 L/min: 0.23 cm H₂O; at 100 L/min: 0.48 cm H₂O
 DEADSPACE VOLUME (mL approx.): **Small** 184 mL **Medium** 193 mL **Large** 209 mL **Extra Large** 220 mL
 Resistance of the anti-asphyxia valve in Single Fault Condition: Inspiratory; 0.8 H₂O per L/s. Expiratory; 0.8 H₂O per L/s
 A-weighted Sound Pressure Level at 1 m distance: 29.98 (A); A-weighted Sound Power Level: 32.99 (A)
 AAV Open to Atmosphere Pressure: 1.45 cm H₂O, AAV Close to Atmosphere Pressure: 1.80 cm H₂O
 Do not expose the mask to temperatures above 140°F (60°C).

BEFORE USE

- Inspect mask daily or prior to each use. Replace the mask if vent holes are blocked, if parts are damaged, or if gel is exposed due to tears or punctures.
- Check the anti-asphyxia valve. With the positive pressure device turned off, verify that the valve flapper is positioned so that room air can flow through the large opening in the valve (Figure 7a). Turn on the system; the flapper should close, and system air should flow through the mask (Figure 7b). If the flapper does not close or does not function properly, replace the mask. Do not block the opening on the anti-asphyxia valve or the vent holes. Ensure that the valve is not blocked with secretions and that the flapper is dry.

VERASEAL® 3 MASK FITTING INSTRUCTIONS

1. Verify the mask is the right size. When using the sizing guide ensure mouth is slightly open. Line the top of the guide to the bridge of the nose and select the smallest size that does not restrict the nose or mouth. The appropriate size may be smaller than expected due to the design.
2. Detach strap by sliding the loop off the Quick Release slot (bottom left strap) (Figure 1), seat the bottom cushion of the mask at the crease of the patient's chin with their mouth slightly open and tip the mask up to contact the patient's nose (Figure 2). Turn the forehead adjustment screw so the forehead pad lightly touches the patient's forehead (Figure 6). The mask may sit lower on the nose than expected, depending on facial features.
3. Slip the headgear over the head, find the detached hook and loop fastener strap and slip it back in the Quick Release slot (Figure 3). Ensure the headgear is flat across the back of the head.
4. Adjust the bottom headgear straps by using the hook and loop fastener tabs (Figure 4).
5. Adjust the upper headgear straps by using the hook and loop fastener tabs: detach the hook and loop fastener tabs and pull straps forward to center mask over face. Ensure headgear is centered by aligning the two straps. Then pull the straps back, adjust strap tension and affix tabs. (Figure 5).
6. Turn the forehead pad adjustment screw to adjust the tension until comfortable. If necessary, readjust top straps for a comfortable fit.
7. If the mask is uncomfortable repeat step 5-7.

NOTE:

1. The mask should not constrict patient's nose or mouth; the patient should be able to breathe freely and comfortably.

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- II. The mask should rest comfortably against the face. If there is excessive pressure on the bridge of the nose, loosen the top straps, slowly turn the adjustment screw counter-clockwise to relieve the pressure on the bridge (Figure 6). Then readjust the top straps and affix the tabs.
8. Connect the mask to the positive pressure device tubing and instruct the patient to breathe normally.
9. If leaks are noticed around the chin or the cheeks, slightly tighten the bottom straps.
10. If leaks are noticed at the bridge of the nose, slowly turn the adjustment screw clockwise to slightly move the pad away from the forehead (Figure 6). Then tighten the top straps slightly until leaks are eliminated.

NOTE:

- I. DO NOT overtighten any straps. Overtightening can worsen leaks. For optimal seal and comfort, mask should be snug, but not tight.

QUICK RELEASE OF MASK AND HEADGEAR

11. Slide hook and loop fastener loop off quick release hook and remove mask assembly.

Veraseal® is a trademark of Sleepnet Corporation. Visit our website at www.sleepnetmasks.com.

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Veraseal® est une marque de commerce de Sleepnet Corporation.