Explanation of Symbols Which Appear on the Device Labeling
Explicación de los símbolos que aparecen en la etiqueta del dispositivo
Erklärung der Symbole auf der Gerätezeichnung
Explication des symboles qui apparaissent sur l'étiquette de l'appareil
Illustrazione dei simboli visualizzati nell'etichetta del dispositivo
Laitteen merkinnöissä näkyvien symbolien selitys
Cihaz Etiketinde Görünen Sembollerin Açıklaması
Explicații cu privire la simbolurile ce apar pe eticheta dispozitivului
Επεξήγηση των συμβόλων που εμφανίζονται στην επισήμανση της συσκευής

Veraseal® 3 Full Face Non Vented Mask

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Veraseal® 3 Full Face Non Vented Mask

CONTENTS: Veraseal® 3 Full Face Non-Vented Mask with headgear
Assembled in USA with US and imported parts
Rx Only

INTENDED USE

The Veraseal® 3 Full Face Non-Vented Mask is a disposable device intended to provide a patient interface for the application of noninvasive ventilation. The mask is to be used as an accessory to ventilators with adequate alarms, and safety systems for ventilator failure, and which are intended to administer positive pressure ventilation. The mask is intended for SHORT TERM SINGLE PATIENT USE (maximum of 7 days) on adult patients (>30kg) who are appropriate candidates for noninvasive ventilation in the hospital or institutional environment.

NOTES

- This mask requires a separate exhalation device.
- The mask system does not contain PVC, DEHP or phthalates.
- This product is not manufactured with natural rubber latex.
- If the patient has ANY reaction to any part of the mask system, discontinue use.
- The mask does not include an anti-asphyxia valve (air entrainment valve).
- Verify that the mask is the right size by using the sizing guide.

CONTRAINDICATIONS

This product should not be used if the patient is experiencing nausea, vomiting, taking a prescription drug that may cause vomiting, or if they are unable to remove the mask on their own.

CAUTIONS

- In the U.S. Federal Law restricts this device for sale by or on the order of a physician.
- Do not expose the mask to temperatures above 140°F (60°C).
- The mask is intended for SINGLE PATIENT USE ONLY. It is intended for short term use (not more than 7 days) and must not be used on multiple persons. Do not sterilize or disinfect.
- Consult attending physician prior to using mask if the patient is using any medications or devices to remove deep secretions.

WARNINGS

- This mask is not suitable for providing life support ventilation.
- This mask must be used with a POSITIVE PRESSURE VENTILATOR WITH AN EXHALATION VALVE. The mask should not be used unless the positive pressure ventilator is turned on and operating properly. When the ventilator is functioning properly, the exhalation valve allows exhaled air to escape to ambient. When the ventilator is not operating, exhaled air may be re-breathed. Re-breathing of exhaled air can, in some circumstances, lead to suffocation.
- Use this mask with baseline pressure of 3 cm H₂O or higher.
- This mask must be used with VENTILATORS WITH EXHALATION VALVES, ADEQUATE ALARMS, AND SAFETY SYSTEMS FOR VENTILATOR FAILURE; use of this mask requires the same level of attention and assistance as the use of a tracheal tube.
- No smoking or open flames may be used when oxygen is in use.
- If oxygen is used with this mask, the oxygen flow must be turned off when the machine is not operating. Explanation of the warning: When the positive pressure device is not in operation, and the oxygen flow is left on, oxygen delivered into the device tubing may accumulate within the positive pressure device enclosure. Oxygen accumulated will create a risk of fire.
- At a fixed flow rate of supplemental oxygen, the inhaled oxygen concentration will vary, depending on pressure settings, patient breathing, and leak rate.
- If the patient has respiratory insufficiency, they should wear this mask ONLY when therapy is being delivered.
- Any unusual skin irritation, chest discomfort, shortness of breath, gastric distension, abdominal pain, belching or flatulence from ingested air, or severe headache experienced by the patient during or immediately after use should be immediately reported to the attending healthcare professional.
- Using a mask may cause tooth, gum or jaw soreness or aggravate an existing dental condition. If symptoms occur, consult a physician or dentist.
- Consult a physician if the patient experiences any of the following symptoms occur while using the mask: drying of the eyes, eye pain, eye infections, or blurred vision. Consult an ophthalmologist if symptoms persist.
- Do not overtighten the headgear straps. Signs of overtightening include: excessive redness, sores, or bulging skin around the edges of the mask. Loosen headgear straps to alleviate symptoms.
- Significant unintentional leak may occur without properly fitting the mask. Follow the MASK FITTING INSTRUCTIONS provided.

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- This mask does not include an anti-asphyxia valve (air entrainment valve) that would allow a patient to breathe if the ventilator failed. This mask is not intended for patients without a spontaneous respiratory drive. This mask should not be used on patients who are uncooperative, obtunded, unresponsive, or unable to remove the mask. This mask may not be suitable for persons with the following conditions: impaired cardiac sphincter function, excessive reflux, impaired cough reflex, and hiatal hernia.

GENERAL INFORMATION

The mask has a service life of 7 days. Please dispose of the mask and accessories in accordance with local and regional guidelines. Per EU MDR 2017/745, please report serious incidents to manufacturer and EU Member State competent authority.

CLEANING AND MAINTENANCE

This mask is disposable and is intended for single patient short term use (maximum of 7 days) only. The mask should be discarded if heavily soiled. The mask may be cleaned daily or as needed using a 70% isopropyl alcohol swab or by wiping with a damp, soapy cloth. If using a damp cloth, use mild soap such as Ivory®. Rinse thoroughly and allow to air dry, out of direct sunlight.

- NOTE:** The gel cushion, although silky to the touch, will not tear, leak, or detach from the mask under normal use. The gel cushion is not indestructible. It will tear, cut, and rip if subjected to mishandling or abuse.

TECHNICAL INFORMATION

Resistance to air flow at 50 L/min: 0.10 cm H₂O Resistance to air flow at 100 L/min: 0.32 cm H₂O Dead Space (approx.): **Small** 180 ml **Medium** 193 ml **Large** 203 ml **Extra Large** 215 ml Operating Temperature: Do not expose the mask to temperatures above 140°F (60°C).

BEFORE USE

- Inspect mask daily or prior to each use. Replace the mask if parts are damaged, or if gel is exposed due to tears or punctures.

VERASEAL® 3 MASK FITTING INSTRUCTIONS

1. Verify the mask is the right size. When using the sizing guide ensure mouth is slightly open. Line the top of the guide to the bridge of the nose and select the smallest size that does not restrict the nose or mouth. The appropriate size may be smaller than expected due to the design.
2. Detach strap by sliding the loop off the Quick Release slot (bottom left strap) (Figure 1), seat the bottom cushion of the mask at the crease of the patient's chin with their mouth slightly open and tip the mask up to contact the patient's nose (Figure 2). Turn the forehead adjustment screw so the forehead pad lightly touches the patient's forehead (Figure 6). The mask may sit lower on the nose than expected, depending on facial features.
3. Slip the headgear over the head, find the detached hook and loop fastener strap and slip it back in the Quick Release slot (Figure 3). Ensure the headgear is flat across the back of the head.
4. Adjust the bottom headgear straps by using the hook and loop fastener tabs (Figure 4).
5. Adjust the upper headgear straps by using the hook and loop fastener tabs: detach the hook and loop fastener tabs and pull straps forward to center mask over face. Ensure headgear is centered by aligning the two straps. Then pull the straps back, adjust strap tension and affix tabs. (Figure 5).
6. Turn the forehead pad adjustment screw to adjust the tension until comfortable. If necessary, readjust top straps for a comfortable fit.
7. If the mask is uncomfortable repeat steps 5-7.

NOTE:

- I. The mask should not constrict patient's nose or mouth; the patient should be able to breathe freely and comfortably.
- II. The mask should rest comfortably against the face. If there is excessive pressure on the bridge of the nose, loosen the top straps, slowly turn the adjustment screw counter-clockwise to relieve the pressure on the bridge (Figure 6). Then readjust the top straps and affix the tabs.
8. Connect the mask to the positive pressure device tubing and instruct the patient to breathe normally.
9. If leaks are noticed around the chin or the cheeks, slightly tighten the bottom straps.
10. If leaks are noticed at the bridge of the nose, slowly turn the adjustment screw clockwise to slightly move the pad away from the forehead (Figure 6). Then tighten the top straps slightly until leaks are eliminated.

NOTE:

- I. DO NOT overtighten any straps. Overtightening can worsen leaks. For optimal seal and comfort, mask should be snug, but not tight.

QUICK RELEASE OF MASK AND HEADGEAR

- Slide hook and loop fastener loop off quick release hook and remove mask assembly.

Veraseal® is a trademark of Sleepnet Corporation.
Visit our website at www.sleepnetmasks.com.

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