


# Veraseal® 3 Full Face Mask

## Quick Start Guide



Available elbow options: Vented, Non-Vented, Non-Vented AAV



Unhook one or both straps and place the mask on the patient's face. Seat the bottom cushion of the mask at the crease of the chin. Tip the mask up to contact the nose, and rethread the bottom strap/s through the slots.



Ensure the headgear is fully extended and lies flat across the back of the head, with the headgear logo facing outward.



Position the upper and lower straps evenly on both sides of the patient's head. Gradually tighten the headgear straps, starting from the bottom and working your way up. Make sure the upper and lower straps are parallel as shown. The mask should fit comfortably snug, but not too tight. Check for any pressure points or discomfort, and adjust as needed.

# Veraseal® 3 Full Face Mask

## Quick Start Guide



Turn the forehead adjustment dial to adjust tension until the patient is comfortable.



If necessary, readjust the straps for a comfortable fit. Check for any pressure points or discomfort, then adjust as needed.

### NOTES:

- Do NOT overtighten the headgear.
- For optimal seal and comfort, the mask should be snug, but not tight.



Verify that the mask is the right size by using the size guide. (Available online)

[https://www.sleepnetmasks.com/wp-content/uploads/2023/11/M90088\\_v0\\_Veraseal3\\_Size\\_Guide.pdf](https://www.sleepnetmasks.com/wp-content/uploads/2023/11/M90088_v0_Veraseal3_Size_Guide.pdf)