Explanation of Symbols
Explicación de símbolos
Erklärung der Symbole
Explication des symboles
Toelichting bij symbolen
Spiegazione dei simboli
Επεξήγηση των συμβόλων

Forklaring på symboler
Förklaring av symboler
Explicação dos símbolos
Wyjaśnienie symboli
Merkkien selitys
Sembollerin açıklaması

Caution or Warning
Precaución o advertencia
Achtung oder Warnung
Avertissement ou mise en garde
Attentiesignaal of waarschuwing
Attenzione o avvertenza
Επισημανση προσοχής ή
προειδοποίηση
Forsigtighed eller advarsel
Försiktighet eller varning
Cuidado ou advertência
Uwaga lub ostrzeżenie
Huomautus tai varoitus
Uyarı ve Önlem

Lot number
Número de lote
Losnummer
Numéro du lot
Partijnummer
Numero di lotto
Αριθμός παρτίδας
Varenummer
Partinummer
Número do lote
Numer partii
Eränumero
Lot numarası

Does not contain natural rubber latex
No contiene látex de caucho natural
Enthält keinen Naturkautschuklatex
Ne contient pas de latex en caoutchouc naturel
Bevat geen latex van natuurlijk rubber
Non contiene lattice di gomma naturale
Δεν περιέχει λάτεξ από φυσικό καουτσούκ
Indeholder ikke naturgummilætex
Innehåller inte rågummilætex
Não contém látex de borracha natural
Nie zawiera lateksu naturalnego.
Ei sisällä luonnonkumista valmistettua lateksia.
Doğal kauçuk lateks içermez.

CE Certification Mark
Marca de certificación CE
CE Zertifizierungsmarke
Marque CE
Marcatura CE
EG-keurmerk
Σήμα πιστοποίησης CE
CE certificeringsmærke
CE-certifieringsmärke
Marca de certificação CE
Oznaczenie certyfikacji CE
CE-sertifiointimerkintä
CE Onay İşareti

Part number
Número de pieza
Teilenummer
Référence
Onderdeelnummer
Numero di articolo
Αριθμός εξαρτήματος
Del nummer
Delnummer
Número da peça
Numer części
Osanumero
Parça numarası

Manufacturer
Fabricante
Hersteller
Fabricant
Fabrikant
Fabbrikante
Κατασκευαστής
Producent
Tillverkare
Fabricante
Producent
Valmistaja
Üretici

Refer to Instructions for Use
Consulte las instrucciones de uso.
Siehe Gebrauchsanleitung
Se reporter au Mode d'emploi
Consultare le Istruzioni per l'uso
Raadpleeg de gebruiksaanwijzing
Ανατρέξτε στις Οδηγίες χρήσης
Der henvises til brugsanvisningen
Se bruksanvisningen
Consulte as instruções de uso
Patrz instrukcje obsługi
Katso käyttöohjeet
Bkz. Kullanma Talimatları

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CONTENTS: **Mojo™** Non Vented Full Face Mask with Headgear
Rx ONLY MADE IN U.S.A.
DOES NOT CONTAIN NATURAL RUBBER LATEX

INTENDED USE

The **Mojo™** Non Vented Full Face Mask is intended to provide a patient interface for application of non-invasive ventilation. The mask is to be used as an accessory to ventilators that have adequate alarms and safety systems for ventilator failure, and which are intended to administer positive pressure ventilation for treatment of respiratory failure or respiratory insufficiency. It is intended for use on adult patients (>30 kg.) who are appropriate candidates for noninvasive ventilation in the hospital or institutional environment.

GENERAL INFORMATION - WARNINGS AND CAUTIONS

CAUTION: U.S. Federal Law restricts this device for sale by or on the order of a physician.

CAUTION: Do not expose the **Mojo™** Non Vented Full Face Mask to temperatures above 122°F (50°C).

CAUTION: This mask may be cleaned and used repeatedly on the same person. Do not sterilize.

WARNING: This mask should be used with a POSITIVE PRESSURE VENTILATOR WITH AN EXHALATION VALVE recommended by the manufacturer or by your physician or respiratory therapist. The mask should not be used unless the positive pressure ventilator is turned on and operating properly. When the ventilator is functioning properly, the exhalation valve allows exhaled air to escape to ambient. When the ventilator is not operating, exhaled air may be re-breathed. Re-breathing of exhaled air can, in some circumstances, lead to suffocation.

WARNING: This mask must be used with ventilators which have adequate alarms and safety systems in the event of a ventilator failure. In the event of ventilator failure, use of this mask requires the same level of attention and assistance as the use of a tracheal tube.

WARNING: No smoking or open flames such as candles may be used when oxygen is in use.

WARNING: If using O₂ enrichment, turn off oxygen supply when not in use. When using a fixed flow rate of supplemental oxygen, the inhaled oxygen percentage will vary, depending on pressure settings, patient breathing pattern, mask size, and leak rate.

WARNING: The **Mojo™** Non Vented Full Face Mask must be used with non-conductive tubing.

WARNING: Do not use the **Mojo™** Non Vented Full Face Mask with baseline pressure less than 3 cm H₂O.

WARNING: To minimize the risk of vomiting during sleep, avoid eating or drinking for three (3) hours before using the mask.

WARNING: This product is indicated for use with POSITIVE PRESSURE VENTILATORS CONTAINING EXHALATION VALVES.

WARNING: This mask does not include an anti-asphyxia valve that would allow a patient to breathe if the ventilator failed. This mask is contraindicated for patients without a spontaneous respiratory drive. This mask should not be used on patients who are uncooperative, obtunded, unresponsive, or unable to remove the mask. This mask may not be suitable for persons with the following conditions: impaired cardiac sphincter function, excessive reflux, impaired cough reflex, and hiatal hernia.

CAUTION: Consult physician prior to using mask if using any medications or devices to remove deep secretions.

WARNING: Any unusual chest discomfort, shortness of breath, gastric distention, abdominal pain, belching or flatulence from ingested air, or severe headache during or immediately after use should be immediately reported to physician.

CAUTION: Any unusual skin irritation should be immediately reported to physician.

Do not use this product if you are experiencing nausea, vomiting, taking a prescription drug that may cause vomiting, or if you are unable to remove the mask on your own.

GENERAL INFORMATION

This mask has a non-transferable three month limited warranty for manufacturer's defects from the date of purchase by the initial consumer. If a mask used under normal conditions should fail, Sleepnet will replace the mask or its components. For additional information about the Sleepnet Warranty, please visit <http://www.sleepnetmasks.com/products/our-promise.cfm>.

NOTE: The mask has a service life of 6 months.

CLEANING AND MAINTENANCE - SINGLE PATIENT USE

Cleaning is recommended on a daily basis. Wash the mask and components in warm water using a mild shampoo or dish detergent. Rinse thoroughly after washing. Visually verify that the swivel elbow is not obstructed and allow the mask to air dry.

NOTE: Inspect mask and valve prior to each use. Replace the mask if parts are damaged or gel is exposed due to tear or punctures.

NOTE: The gel cushion, although silky to touch, will not tear, leak or detach from the mask under normal use. Follow recommended cleaning and maintenance protocols. The gel cushion is not indestructible. It will tear, cut and rip if subjected to mishandling or abuse. Do not clean in a dishwashing machine. The warranty does not cover these types of product failures. Therefore, please handle your new **Mojo™** Non Vented Full Face Mask with care.

TECHNICAL INFORMATION

DEADSPACE VOLUME (mL approx.) **Small** 165 mL **Medium** 190 mL **Large** 210 mL
Extra Large 215 mL

USER INSTRUCTIONS FOR **MOJO™** HEADGEAR

1. Unhook the right or left magnetic clip from the mask.
2. Place the mask and headgear assembly on head with blue fabric facing outward. Hook the magnetic clip in place.
3. Adjust upper and lower headgear straps; unhook the Velcro® tabs and pull forward to center mask over face, and adjust strap tension and length. Affix tabs.

NOTE: The mask should rest comfortably against the face and still provide an effective air seal. Straps should be adjusted to stabilize mask; pulling the straps too tight can cause skin irritation. The mask should not constrict your nose or mouth; you should be able to breathe freely and comfortably. If the mask and headgear are adjusted correctly and a problem persists, contact your healthcare professional.

4. Attach swivel adapter to device tubing (bi-level, ventilator); turn device on and verify airflow through mask.
5. Gently bend pliable ring embedded inside mask shell to change mask shape, create a comfortable fit, and reduce leaks.
6. Rotate the tensioning screw to eliminate leaks at the bridge of the nose. Do not attempt to decrease bridge leaks by over-tightening the upper strap.

QUICK RELEASE OF MASK AND HEADGEAR

7. Grasp lower strap and lift away from face. Remove mask and headgear assembly.

Mojo™ is a trademark of Sleepnet Corporation.

Visit our website at www.sleepnetmasks.com

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SOLAMENTE PRESCRIZIONE MEDICA PRODOTTO IN U.S.A.