

SLEEPNET MASK PRESSURE FLOW CHART

Pressure (cm of H ₂ O)	Approx. Exhaust flow (lpm)			
	IQ Nasal mask	Phantom Nasal mask	Minime Pediatric mask	Mojo Full Face mask
3	12.4	9.5	14.0	21.2
4	14.2	11.9	16.1	24.1
5	15.5	13.2	18.3	26.8
6	17.3	15.0	19.8	29.0
7	19.0	16.2	20.4	31.0
8	20.3	17.4	23.2	33.3
9	21.8	18.7	24.3	35.3
10	23.0	19.6	25.6	37.2
11	24.1	20.8	26.8	38.9
12	25.0	21.8	27.9	41.1
13	25.8	22.9	28.8	42.5
14	27.9	23.9	29.8	44.5
15	29.5	24.6	30.7	45.9
16	30.0	25.1	31.9	47.5
17	31.0	25.9	32.8	48.8
18	31.5	26.8	33.7	50.5
19	31.9	27.5	34.3	51.8
20	32.5	28.5	35.0	53.1