

Sleepnet®

Mojo®

featuring breakthrough AIRgel™ technology

## Full-Face Mask Size Gauge


SMALL 50835

MEDIUM 50836

LARGE 50837

To measure, hold size gauge up to patient's face and measure from indent on nose bridge to indent under lips. Relax jaw with mouth slack and slightly open (about 1" / 2.5cm).

*Example:* Patient's size is SMALL



For ordering information call **1.800.742.3646**  
or visit [www.sleepnetmasks.com](http://www.sleepnetmasks.com) to learn more.

**Sleepnet Corporation**

5 Merrill Industrial Drive, Hampton, New Hampshire 03842 USA