

Our MiniMe® Vented and Non-Vented masks feature a patented AIR°gel cushion and flexible shell with Custom Fit Technology™. MiniMe is designed for children and FDA cleared for ages 2 – 12. As the only customizable mask on the market, MiniMe provides your smallest patients with uncompromising fit and maximum comfort.

## The Smaller Side of Comfort

### Custom Fit Technology

The shapeable shell of each mask can be easily molded to fit individual needs. Flex, fit, forget.



### Clearest Field of Vision Possible

Minimal visual interference supports therapy compliance. Children can watch TV, wear glasses, or read a book.

### AIR°gel

Maximum comfort, minimum weight. Our patented AIR°gel minimizes pressure points and skin irritation.

### Vented and Non-Vented Versions Available

For CPAP and Bi-level therapy in children.

### Quick-Release Headgear

Lightweight, breathable headgear slides easily on and off so kids can manage this task on their own.

### One Convenient Design

Engineered to fit effectively for pediatric patients ages 2 – 12.

### Advanced Cushion Technology (A.C.T.)

Our ergonomic cushion technology delivers an effective seal for the most comfortable fit.



Made in the USA

# What makes our MiniMe Mask different?

- FDA Cleared for ages 2 – 12. As the leading mask designed just for children, MiniMe has a historical track record of success in the pediatric market.
- Custom Fit Technology is unique to Sleepnet. Our soft-shell mask has a shapeable wire allowing for optimal flexibility, fit and comfort.
- One-size design with quality construction eliminates the need for replacement parts and reduces SKUs.

## About Advanced Cushion Technology (A.C.T.)

Our patented AIR<sup>o</sup>gel with an ergonomic streamlined design and luxurious feel reduces pressure while maintaining an effective seal.



## Custom Fit Technology

BEND IT.



MOLD IT .




SHAPE IT .



## MiniMe Mask Specifications:

Product Number	Mask Style	Leak Port	Headgear Size	Mask Size	Mask Weight	Mask Dimensions
60213	Vented	✓	Small	One size	31.5 g	1.60"D x 3.10"W x 2.14"H
60214	Vented	✓	Medium	One size	31.5 g	1.60"D x 3.10"W x 2.14"H
60215	Vented	✓	Large	One size	31.5 g	1.60"D x 3.10"W x 2.14"H
60253	Non-Vented	✗	Small	One size	31.5 g	1.59"D x 3.10"W x 2.18"H
60254	Non-Vented	✗	Medium	One size	31.5 g	1.59"D x 3.10"W x 2.18"H
60255	Non-Vented	✗	Large	One size	31.5 g	1.59"D x 3.10"W x 2.18"H



 Does not contain natural rubber latex

# Aura<sup>®</sup>

Your Mask, Your Way



## Contact us today.

To place your order, call an authorized distributor.

Call 1.800.742.3646

Visit [www.sleepnetmasks.com](http://www.sleepnetmasks.com) to learn more.

Sleepnet Corporation, 5 Merrill Industrial Drive, Hampton, New Hampshire 03842 USA

