


MiniMe Headgear Sizing Chart

Head Size - Metric	35.5	38	41	43	46	48.5	49.5	52
Head Size - Inches	14	15	16	17	18	19	19 1/2	20 1/4
60215 Large Headgear only # 55243					○	●	●	●
60214 Medium Headgear only # 55242			○	●	●	○		
60213 Small Headgear only # 55241	○	●	●	○				



To determine head size, measure circumference of head, just above eyebrows as if for a sweat band. Select the most appropriate size.

● Recommended
○ May also fit


Breathe different.

- Patented AIR[°]gel cushion delivers an effective seal with no irritating pressure points on the skin.
- FDA Cleared for ages 2 – 12. As the leading mask designed just for children, MiniMe has a historical track record of success in the pediatric market.
- Custom Fit Technology™ is unique to Sleepnet®. Our soft-shell mask has a shapeable wire allowing for optimal flexibility, fit and comfort.
- Lightweight, breathable headgear slides easily on and off so kids can manage this task on their own.
- One-size design with quality construction eliminates the need for replacement parts and reduces SKUs.



CE 0123

Does not contain natural rubber latex



Made in the USA

Call 1.800.742.3646 or visit www.sleepnetmasks.com to learn more.

Sleepnet Corporation, 5 Merrill Industrial Drive, Hampton, New Hampshire 03842 USA

Scan code to view MiniMe Fitting Tips video.

