

The Mojo 2 Vented and Non-Vented masks feature a patented AIR<sup>o</sup>gel cushion, moldable shell with Custom Fit Technology™ and a Forehead Adjustment Screw. Together, they provide a personalized fit with maximum comfort.

## AIR<sup>o</sup>gel<sup>®</sup>

Maximum comfort, minimum weight. Our patented AIR<sup>o</sup>gel minimizes pressure points and skin irritation.



## Custom Fit Technology

The shapeable shell of each mask can be easily molded to fit individual needs.

**Bend It. Mold It. Shape It.**

## Oxygen Port

Allows easy delivery of supplemental oxygen.

## 360° Rotating Elbow

For optimal tube positioning and easier mobility.

## Forehead Adjustment Screw

Removes pressure from the bridge of the nose for personalized fit.

## AIR<sup>o</sup>gel<sup>®</sup> Forehead Cushion

## Vented, Non-Vented, and Non-Vented AAV Versions Available

## Quick-Release Magnetic Headgear Clips

Lightweight, breathable headgear disconnects from mask easily with one hand.

## Size Indicator

Available in four sizes

**S M L XL**



Made in the USA  
of US and imported parts

# What makes our Mojo 2 Mask different?

- Custom Fit Technology - a unique soft-shell mask with a shapeable wire allowing for optimal flexibility, fit, and comfort.
- Adjustment Screw - provides a wide, continuous range of adjustment ensuring a custom fit while removing pressure from the forehead.
- Available in four sizes with magnetic connectors and simple, breathable quick-release headgear.

## About AIR°gel

Patented AIR°gel cushion delivers an effective seal while minimizing pressure points on the skin.



## Custom Fit Technology

BEND IT.



MOLD IT.



SHAPE IT.



## Mojo 2 Masks with Headgear Specifications:

Product Number	Style	Size	Leak Port	Entrianment Valve	Mask Weight	Mask Dimensions
50817	Vented	Small	✓	✓	133.8g	3.60"D x 4.30"W x 7.25"H
50818	Vented	Medium	✓	✓	141g	3.60"D x 4.40"W x 7.40"H
50819	Vented	Large	✓	✓	144.5g	3.60"D x 4.50"W x 7.60"H
50820	Vented	Extra Large	✓	✓	146.5 g	3.60"D x 4.55"W x 7.80"H
50912	Non-Vented	Small	✗	✗	133.2g	3.70"D x 4.30"W x 7.25"H
50913	Non-Vented	Medium	✗	✗	140.4g	3.70"D x 4.40"W x 7.40"H
50914	Non-Vented	Large	✗	✗	143.9g	3.70"D x 4.50"W x 7.60"H
50915	Non-Vented	Extra Large	✗	✗	145.9 g	3.70"D x 4.55"W x 7.80"H
50941	AAV Non-Vented	Small	✗	✓	133.7g	3.65"D x 4.30"W x 7.25"H
50942	AAV Non-Vented	Medium	✗	✓	140.9g	3.65"D x 4.40"W x 7.40"H
50943	AAV Non-Vented	Large	✗	✓	144.4g	3.65"D x 4.50"W x 7.60"H
50944	AAV Non-Vented	Extra Large	✗	✓	146.4g	3.65"D x 4.55"W x 7.80"H



Does not contain natural rubber latex

## Contact us today.

General inquiries email  
[info@sleepnetcorp.com](mailto:info@sleepnetcorp.com) or  
 call 1 (603) 758-6600

Visit [www.sleepnetmasks.com](http://www.sleepnetmasks.com) to learn more.

Sleepnet Corporation, 5 Merrill Industrial Drive,  
 Hampton, New Hampshire 03842 USA