




# iQ<sup>®</sup> Blue

featuring AIR<sup>°</sup>gel<sup>®</sup> with Advanced Cushion Technology

Introducing the next generation in Airgel. Our iQ<sup>®</sup> Blue Nasal Vented and Non-Vented masks feature our patented AIR<sup>°</sup>gel with *Advanced Cushion Technology* and flexible shell with *Custom Fit Technology*<sup>™</sup>. Together, they provide an optimal fit with maximum comfort. Unlike masks with rigid shells, iQ Blue adjusts to truly fit the contours of your patients' face.

## The Smart Choice for a Custom Fit



# What makes our iQ Blue Mask different?

- Custom Fit Technology - a unique soft-shell mask with a shapeable wire allowing for optimal flexibility, fit, and comfort.
- Our AIR°gel with Advanced Cushion Technology creates an ergonomic gel cushion that reduces pressure points on the skin.
- Clearest field of vision possible - users can even wear glasses with this mask.
- No forehead pad, yet retains outstanding stability.

# About Advanced Cushion Technology(A.C.T.)

Our patented AIR°gel with an ergonomic design and luxurious feel reduces pressure while maintaining an effective seal.



## Custom Fit Technology

BEND IT.



MOLD IT.



SHAPE IT.



## iQ Blue Mask Specifications:

Product Number	Mask Style	Leak Port	Headgear Style	Mask Size	Tube Length	Mask Weight	Mask Dimensions
50647	Vented	✓	3-Point	OS	13"	63.8 g	3.04"D x 2.95"W x 2.60"H
50655	Vented	✓	StableFit	OS	13"	63.8 g	3.04"D x 2.95"W x 2.60"H
50250	Vented	✓	3-Point	OS	6"	63.8 g	3.04"D x 2.95"W x 2.60"H
50255	Vented	✓	StableFit	OS	6"	63.8 g	3.04"D x 2.95"W x 2.60"H
50660	Non-Vented	✗	3-Point	OS	6"	63.9 g	3.08"D x 3.13"W x 2.50"H

CE 0123



Does not contain natural rubber latex

# Aura®

Your Mask, Your Way



## Contact us today.

To place your order, call an authorized distributor.

Call 1.800.742.3646  
Visit [www.sleepnetmasks.com](http://www.sleepnetmasks.com) to learn more.  
Sleepnet Corporation, 5 Merrill Industrial Drive,  
Hampton, New Hampshire 03842 USA



Scan code to view iQ Fitting Tips video.