



# Innova<sup>®</sup> Nasal

featuring AIR<sup>°</sup>gel<sup>®</sup> with Advanced Cushion Technology

The Innova<sup>®</sup> Nasal CPAP mask features our patented AIR<sup>°</sup>gel with Advanced Cushion Technology and Active Headgear Connector<sup>™</sup>. Together, they provide an optimal fit, maximum comfort, and freedom to move without breaking the seal.

## Move Freely with Our New Nasal Mask

### AIR<sup>°</sup>gel

Maximum comfort, minimum weight. Our patented AIR<sup>°</sup>gel minimizes pressure points and skin irritation.



### Touchless Spacebar<sup>™</sup>

Unique design eliminates any contact with the patient's forehead.

### Size Indicator

available in two sizes.

**S/M** **M/L**

### Advanced Cushion Technology (A.C.T.)

Our ergonomic cushion technology delivers an effective seal for the most comfortable fit.

### Quick-Release Headgear

Lightweight, breathable headgear slides off easily with one hand.

### Active Headgear Connector

Moves as you move for a reliable seal.

### 360° Rotating Elbow

For optimal tube positioning and easier mobility



Made in the USA

## What makes our Innova Mask different?

- Our unique *Active Headgear Connector* moves when you move to help maintain a seal throughout the night.
- Easy-to-use with a *Quick-Release Headgear Connector*.
- *Touchless Spacebar™* eliminates any contact with the forehead, yet retains stability.
- Available in two sizes with 4-point simple and breathable headgear.

### About AIR°gel

Patented AIR°gel cushion delivers an effective seal while minimizing pressure points on the skin.



## Active Headgear Connector



Fifteen degrees of rotation provide relief from headgear tension, allowing the mask to move with the patient for a reliable seal.

### Innova Mask with Headgear Specifications:

Product Number	Style	Leak Port	Size	Mask Weight	Mask Dimensions
50162	Vented	✓	S/M	57.0 g	3.39"D x 2.66"W x 2.27"H
50163	Vented	✓	M/L	65.1 g	3.77"D x 2.99"W x 2.67"H
50592	Non-Vented	✗	S/M	57.0 g	3.39"D x 2.66"W x 2.27"H
50593	Non-Vented	✗	M/L	65.1 g	3.77"D x 2.99"W x 2.67"H
55051	Headgear	N/A	OS	N/A	N/A

 Does not contain natural rubber latex

 0123

### Also from Sleepnet:

## Mojo®

Your Mask, Your Way

Our full face mask also features AIR°gel.



### Contact us today.

To place your order, call an authorized distributor.

Call 1.800.742.3646

Visit [www.sleepnetmasks.com](http://www.sleepnetmasks.com) to learn more.

Sleepnet Corporation, 5 Merrill Industrial Drive, Hampton, New Hampshire 03842 USA

