

Collaborating with leading children's hospitals, Sleepnet now offers the MiniMe 2 mask in **4** sizes and **new** soft, 5-adjustment point headgear in **2** sizes. FDA cleared for children ages 2 - 12, Sleepnet's MiniMe 2 features a proprietary shapeable pediatric mask with our patented AIR°gel technology. The **new** MiniMe 2 provides your smallest patients the ultimate fit and comfort.



The Next Generation of Pediatric Masks

AIR°gel®

Maximum comfort, minimum weight. Our patented AIR°gel makes the MiniMe 2 the most gentle fitting mask on the market.



Touchless Spacebar™

Unique design eliminates any contact with the patient's forehead.

Size Indicator

Both Vented and Non-Vented come in four sizes: **XS, S, M, L**

Custom Fit Technology

The shapeable shell mask can be easily molded to fit individual needs.
Flex, Fit, Forget.

Advanced Cushion Technology (A.C.T.)

Our ergonomic cushion technology delivers an effective seal for the most comfortable fit.

360° Rotating Elbow

For optimal tube positioning and easier mobility

Vented and Non-Vented Versions Available

CPAP, BPAP and Ventilation Therapy in children



Made in the USA
of US and imported parts

What makes our MiniMe 2 Mask different?

- The MiniMe 2 incorporates our softest and lightest gel technology providing superior comfort and seal for all.
- Custom Fit Technology - a unique soft shell shapeable mask allowing for optimal flexibility, fit, and comfort.
- Mask available in **4** sizes.
- New simple and breathable headgear with 5-adjustment points available in **2** sizes.

Advanced Cushion Technology (A.C.T.)


Our patented AIR^ogel with an ergonomic design and luxurious feel reduces pressure while maintaining an effective seal.



Custom Fit Technology



Flex



Fit



Forget

MiniMe 2 Masks with Headgear Specifications:

Product Number	Style	Headgear size	Mask Size	Mask Weight	Dead Space	Mask Dimensions
60216	Vented	XSmall/Small	XS	33.5 g	20 ml	2.65"D x 3.12"W x 3.94"H
60217	Vented	XSmall/Small	S	31.6 g	25 ml	2.65"D x 3.12"W x 3.94"H
60218	Vented	Medium/Large	M	43.2 g	37 ml	2.81"D x 3.12"W x 4.16"H
60221	Vented	Medium/Large	L	40.3 g	40 ml	2.81"D x 3.12"W x 4.16"H
60256	Non-Vented	XSmall/Small	XS	33.5 g	55 ml	2.65"D x 3.12"W x 3.94"H
60257	Non-Vented	XSmall/Small	S	31.6 g	60 ml	2.65"D x 3.12"W x 3.94"H
60258	Non-Vented	Medium/Large	M	43.2 g	72 ml	2.81"D x 3.12"W x 4.16"H
60261	Non-Vented	Medium/Large	L	40.3 g	75 ml	2.81"D x 3.12"W x 4.16"H
55055	Headgear	XSmall/Small	N/A	N/A	N/A	N/A
55056	Headgear	Medium/Large	N/A	N/A	N/A	N/A

MiniMe 2 Headgear Sizing Chart:

Head Size - Metric		35.5	38	41	43	46	48	51	53.5	56	58.5	61	63.5
Head Size - Inches		14	15	16	17	18	19	20	21	22	23	24	25
55055	XSmall/Small Headgear												
55056	Medium/Large Headgear												

To determine head size, measure circumference of head, just above eyebrows as if for a sweat band. Select the most appropriate size.

- Recommended fit

- Fits in some cases

Two generations of pediatric masks:

MiniMe[®]



MiniMe[®] 2



Contact us today.

To place your order, call an authorized distributor.

Call 1.800.742.3646

Visit www.sleepnetmasks.com to learn more.

Sleepnet Corporation, 5 Merrill Industrial Drive, Hampton, New Hampshire 03842 USA

