



Innova[®] Hospital


featuring AIR[°]gel[®] with Advanced Cushion Technology

The Innova[®] Hospital disposable mask features our patented AIR[°]gel with Advanced Cushion Technology and Active Headgear Connector[™]. Together, they provide an optimal fit, maximum comfort, and freedom to move without breaking the seal.

Move Freely with Our Disposable Hospital Nasal Mask

AIR[°]gel

Maximum comfort, minimum weight. Our patented AIR[°]gel minimizes pressure points and skin irritation.



Touchless Spacebar[™]

Unique design eliminates any contact with the patient's forehead.

Size Indicator

available in two sizes.

S/M **M/L**

Advanced Cushion Technology (A.C.T.)

Our ergonomic cushion technology delivers an effective seal for the most comfortable fit.

Quick-Release Headgear

Lightweight, breathable headgear slides off easily with one hand.

Active Headgear Connector


Moves as you move for a reliable seal.

Vented and Non-Vented Versions Available

For CPAP, Bi-level, and ventilator therapy

360° Rotating Elbow

For optimal tube positioning and easier mobility



Made in the USA

100% designed and manufactured in the USA • Made with environmentally conscious production practices • Our products are backed by our iron-clad warranty against manufacturing defects.

What makes our Innova Mask different?

- Vented and nonvented versions are available for easy hook up to CPAP, Bi-level, or ventilaor therapy..
- Easy-to-use with a *Quick-Release Headgear Connector*.
- *Touchless Spacebar™* eliminates any contact with the forehead, yet retains stability.
- Available in two sizes with 4-point simple and breathable headgear.

Advanced Cushion Technology(A.C.T.)

Our patented AIR°gel with an ergonomic design and luxurious feel reduces pressure while maintaining an effective seal.



Active Headgear Connector



Fifteen degrees of rotation provide relief from headgear tension, allowing the mask to move with the patient for a reliable seal.

Innova Hospital Mask with Headgear Specifications:

Product Number	Style	Leak Port	Size	Mask Weight	Mask Dimensions
50180	Vented	✓	S/M	57.0 g	3.39"D x 2.66"W x 2.27"H
50181	Vented	✓	M/L	65.1 g	3.77"D x 2.99"W x 2.67"H
50594	Non-Vented	✗	S/M	57.0 g	3.39"D x 2.66"W x 2.27"H
50595	Non-Vented	✗	M/L	65.1 g	3.77"D x 2.99"W x 2.67"H



Does not contain natural rubber latex

Also from Sleepnet:

Veraseal 2

A Giant Step Forward in Patient Comfortt



Contact us today.

To place your order, call an authorized distributor.

Call 1.800.742.3646

Visit www.sleepnetmasks.com to learn more.

Sleepnet Corporation, 5 Merrill Industrial Drive, Hampton, New Hampshire 03842 USA

