

The V3 combines AIR°gel® technology with acute care performance for a revolutionary home care experience. The V3 fits quickly and easily while enhancing patient comfort.

## AIR°gel®

Maximum comfort, minimum weight. Our patented AIR°gel minimizes pressure points and skin irritation.

## AIR°gel® Forehead Cushion

## Forehead Adjustment Screw

## Quick-Release Slot

Lightweight, breathable headgear slides off easily with one hand.

## 360° Rotating Elbow

For optimal tube positioning and easier mobility.

## Oxygen Port

Allows easy delivery of supplemental oxygen.

## Size Indicator

Available in four sizes

**S M L XL**



# What makes our V3 Mask different?

- Adjustment Screw -provides a wide, continuous range of adjustment ensuring a custom fit while removing pressure from the bridge of the nose.
- AIR°gel Forehead Pad - stabilizes the mask comfortably on the patient's face during a night's sleep minimizing pressure points.
- Available in four sizes with simple and breathable quick-release headgear.

## Featuring Patented AIR°gel<sup>®</sup> Technology

Our patented AIR°gel with an ergonomic design and luxurious feel reduces pressure while maintaining an effective seal.



## V3 Masks with EZfit Headgear For Comfortable & Secure Treatment at Night



### V3 Product Specifications:

Product Number	Style	Leak Port	Entrianment Valve	Size	Mask Weight	Mask Dimensions
50863	Vented	√	√	Small	109.2g	3.60"D x 4.10"W x 7.25"H
50864	Vented	√	√	Medium	114.3g	3.60"D x 4.30"W x 7.40"H
50865	Vented	√	√	Large	119.4g	3.60"D x 4.38"W x 7.60"H
50866	Vented	√	√	Extra Large	120.1 g	3.60"D x 4.42"W x 7.80"H



Does not contain natural rubber latex

Contact us today.

General inquiries email  
[info@sleepnetcorp.com](mailto:info@sleepnetcorp.com) or call 1 (603) 758-6600

Visit [www.sleepnetmasks.com](http://www.sleepnetmasks.com) to learn more.

Sleepnet Corporation, 5 Merrill Industrial Drive,  
 Hampton, New Hampshire 03842 USA