

Mojo® 2

Your mask, your way

The Mojo® 2 provides a personalized fit with maximum comfort. Our patented AIR°gel cushion and Custom Fit Technology allows the user to mold the Mojo® 2 for an individualized fit.

Flex, Fit, Forget

Forehead Adjustment Dial

Removes pressure from the bridge of the nose which is an essential feature for comfort and compliance.



AIR°gel®
Forehead Cushion

Custom Fit Technology

Patented technology allows each user the ability to accommodate for various facial structure, changes in weight, and facial asymmetry.

Sizes Available

S M L XL

AIR°gel®

Maximum comfort and minimal headgear tension. AIR°gel® minimizes pressure points and skin irritation.

Oxygen Port

Allows easy delivery of supplemental oxygen.

Magnetic Headgear Clips

Easy on and off.

360° Rotating Elbow

For optimal tube positioning and easier mobility.

Available elbow options

Vented
Non-Vented
Non-Vented AAV

What makes our Mojo® 2 Mask different?

- Custom Fit Technology** - A unique soft-shell mask which allows for optimal flexibility, fit, and comfort.
- Adjustment Dial** - Provides a wide, continuous range of adjustment ensuring a custom fit while removing pressure from the nasal bridge.
- Magnetic Headgear Clips** - Easy on and off.
- Available in four sizes** - Vented, Non-Vented, and Non-Vented AAV.

S M L XL



AIRgel® Cushion Technology

Our patented AIRgel technology with an ergonomic design and luxurious feel reduces pressure while maintaining an effective seal.

Color coded elbow options for easy identification to meet your homecare and hospital needs.



Vented Mask with an anti-asphyxia valve for use with CPAP or bi-level devices.



Non-Vented AAV Mask with an anti-asphyxia valve for use with single limb circuits.



Non-Vented Mask for use with dual limb circuits.

Mojo® 2 Masks with Headgear Specifications:

Product Number	Description	Size	Leak Port	Entrainment Valve
50817	Mojo® 2 Full Face Moldable Mask Vented with Headgear	Small	✓	✓
50818	Mojo® 2 Full Face Moldable Mask Vented with Headgear	Medium	✓	✓
50819	Mojo® 2 Full Face Moldable Mask Vented with Headgear	Large	✓	✓
50820	Mojo® 2 Full Face Moldable Mask Vented with Headgear	Extra Large	✓	✓
50912	Mojo® 2 Full Face Moldable Mask Non-Vented with Headgear	Small	X	X
50913	Mojo® 2 Full Face Moldable Mask Non-Vented with Headgear	Medium	X	X
50914	Mojo® 2 Full Face Moldable Mask Non-Vented with Headgear	Large	X	X
50915	Mojo® 2 Full Face Moldable Mask Non-Vented with Headgear	Extra Large	X	X
50941	Mojo® 2 Full Face Moldable Mask AAV Non-Vented with Headgear	Small	X	✓
50942	Mojo® 2 Full Face Moldable Mask AAV Non-Vented with Headgear	Medium	X	✓
50943	Mojo® 2 Full Face Moldable Mask AAV Non-Vented with Headgear	Large	X	✓
50944	Mojo® 2 Full Face Moldable Mask AAV Non-Vented with Headgear	Extra Large	X	✓
55051	EZFit Headgear	One Size		

Also by Sleepnet®
Ascend® Full Face
Replaceable AIRgel® cushion mask



Contact and support
Sleepnet representative in your area
www.sleepnetmasks.com
Phone: 1-800-742-3646

The following product is protected by
US patents 10,974,007, 9,199,054, 7,938,117,
D892,304 and other patents pending



Does not contain
natural rubber latex



sleepnetmasks.com