

MiniMe® 2

The smaller side of comfort

Collaborating with leading children's hospitals. Sleepnet's MiniMe® 2 pediatric nasal mask offers a softer approach to therapy. The MiniMe® 2 offers a new level of flexibility by using exclusive Custom Fit Technology. Sleepnet's patented AIR°gel cushion provides pediatric patients the ultimate fit, comfort, and seal. FDA cleared for children 2-12.



Custom Fit Technology

The shapeable shell mask can be easily molded for a personal fit.

Flex, Fit, Forget!

AIR°gel®


Maximum comfort and minimal headgear tension. AIR°gel® minimizes pressure points and skin irritation.

Touchless Spacebar

Unique design provides superior stability, holding a seal without contacting the skin.

Available sizes:

XS S M L



360° Rotating Elbow

For optimal tube positioning and easier mobility.

Vented and Non-Vented Versions Available

MiniMe® 2 can be used for both hospital and homecare.

What makes our MiniMe® 2 mask different?

- Custom Fit Technology** - Mold the soft shell mask for a personalized fit.
- Touchless Spacebar** - Unique design provides superior stability, holding a seal without contacting the skin.
- Mask available in 4 sizes
 - XS
 - S
 - M
 - L
- Headgear** - Breathable 5-point headgear, available in two sizes
- FDA cleared for ages 2-12



AIRgel® Cushion Technology

Our patented AIRgel technology with an ergonomic design and luxurious feel reduces pressure while maintaining an effective seal.

Custom Fit Technology



FLEX



FIT



FORGET

MiniMe® 2 masks with headgear specifications:

Product Number	Description	Headgear Size	Mask Size	Dead Space
60216	MiniMe® 2 Pediatric Vented Nasal Mask with Headgear	XSmall / Small	XSmall	20 ml
60217	MiniMe® 2 Pediatric Vented Nasal Mask with Headgear	XSmall / Small	Small	50 ml
60218	MiniMe® 2 Pediatric Vented Nasal Mask with Headgear	Medium / Large	Medium	37 ml
60221	MiniMe® 2 Pediatric Vented Nasal Mask with Headgear	Medium / Large	Large	40 ml
60256	MiniMe® 2 Pediatric Non-Vented Nasal Mask	XSmall / Small	XSmall	55 ml
60257	MiniMe® 2 Pediatric Non-Vented Nasal Mask	XSmall / Small	Small	60 ml
60258	MiniMe® 2 Pediatric Non-Vented Nasal Mask	Medium / Large	Medium	72 ml
60261	MiniMe® 2 Pediatric Non-Vented Nasal Mask	Medium / Large	Large	75 ml
55055	MiniMe® 2 Headgear only	XSmall / Small	N/A	N/A
55056	MiniMe® 2 Headgear only	Medium / Large	N/A	N/A

MiniMe® 2 headgear sizing chart:

	■ - Recommended Fit								■ - Fits in some cases			
Head Size - Metric	35.5	38	41	43	46	48	51	53.5	56	58.5	61	63.5
Head Size - Inches	14	15	16	17	18	19	20	21	22	23	24	25
55055	■	■	■	■	■	■	■	■				
55056						■	■	■	■	■	■	■

Also by Sleepnet®

Phantom® 2

featuring breakthrough AIRgel® technology



Contact and support

Sleepnet representative in your area
www.sleepnetmasks.com
Phone: 1-800-742-3646

The following product is protected by
US patents 10,974,007, 9,199,054, 7,938,117,
D892,304 and other patents pending



Does not contain
natural rubber latex



sleepnetmasks.com