

# Aura<sup>®</sup>

## Enter a new aura of comfort and seal

Enjoy comfort and a great seal with the world's first nasal mask featuring both Custom Fit Technology and an Active Headgear Connector. With three sizes and a Touchless Spacebar, Aura will give you the rest you deserve.



### Custom Fit Technology

Patented technology allows each user the ability to accommodate for various facial structure, changes in weight, and facial asymmetry.

### AIR<sup>®</sup>gel

Maximum comfort and minimal headgear tension. AIR<sup>®</sup>gel minimizes pressure points and skin irritation.

### Quick Release Tab

Lightweight breathable headgear slides off easily with one hand.

### 360° Rotating Elbow

For optimal tube positioning and easier mobility.

### Touchless Spacebar

Unique bendable design provides superior stability, holding a seal without contacting the skin.

### Size Indicator

Available in three sizes.



### Advanced Cushion Technology (A.C.T.)

Our ergonomic cushion technology delivers an effective seal for the most comfortable fit.

### Active Headgear Connector

Moves as patient moves for a reliable seal.



## What makes our Aura Mask different?

**AIR°gel® Cushion Technology** - Maximum comfort and minimal headgear tension. AIR°gel® minimizes pressure points and skin irritation.

**Custom Fit Technology** - Patented technology allows each user the ability to accommodate for various facial structure, changes in weight, and facial asymmetry.

**Touchless Spacebar** - Unique bendable design provides superior stability, holding a seal without contacting the skin.

**Three sizes available** - Simple and breathable quick release EZFit Headgear.



### AIR°gel® Cushion Technology

Our patented AIR°gel technology with an ergonomic design and luxurious feel reduces pressure while maintaining an effective seal.

## Custom Fit Technology



**FLEX**



**FIT**



**FORGET**

### Aura Mask and Headgear Specifications:

Product Number	Case Quantity	Style	Leak Port	Size	Mask Weight	Mask Dimensions
50176	50	Vented	✓	S	61.5g	3.80"D x 2.75"W x 5.50"H
50177	50	Vented	✓	M	62.5g	3.80"D x 2.85"W x 5.50"H
50178	50	Vented	✓	L	71.5 g	3.90"D x 3.20"W x 5.50"H
55051	10	EZFit Headgear	N/A	One Size	N/A	N/A

Also by Sleepnet®

**Mojo®**

featuring breakthrough AIR°gel® technology



### Contact us today:

To place your order, call an authorized distributor.

Call 1-800-742-3646

Visit [www.sleepnetmasks.com](http://www.sleepnetmasks.com) to learn more

The following product is protected by US patents 10,974,007, 9,199,054 and other patents pending



Does not contain natural rubber latex



[sleepnetmasks.com](http://sleepnetmasks.com)