



**PEDIATRIC**

# MiniMe<sup>®</sup> Headgear Sizing Chart

Headgear Size - Metric	35.5	38	41	43	46	48.5	49.5	52
Head Size - Inches	14	15	16	17	18	19	19 1/2	20 1/2
60215 Large Headgear only #55243					○	●	●	●
60214 Medium Headgear only #55242			○	●	●	○		
60213 Small Headgear only #55241	○	●	●	○				



Recommended - ●      May also fit - ○



To determine head size, measure circumference of head, just above eyebrows as if for a sweat band. Select the most appropriate size.

## Big comfort for small patients

**Comfort** - Sleepnet's AIR<sup>o</sup>gel<sup>®</sup> is a leading edge innovation developed and designed to enhance the comfort and effectiveness of our respiratory masks.

**Personalized fit** - Custom Fit Technology allows the user to easily shape the mask by gently bending and molding it along the sides and top to provide a consistent and secure fit.

**Easy to fit** - 3 Headgear options available: (S,M,L); One size cushion design.

**FDA Cleared for ages 2-12** - Specifically engineered with children in mind, Sleepnet's MiniMe<sup>®</sup> provides unmatched comfort, a non-invasive design for easy adjustment, and a pediatric size frame for the extra small faces we care so much about.



Does not contain natural rubber latex

