




Innova® Full Face

featuring AIR^ogel[®] with Advanced Cushion Technology

It's time for a better night's sleep. Our lightweight full face mask has our Active Headgear Connector™ that moves when you move to prevent leaks. The sleek design of the Flexible Spacebar™ eliminates contact with your forehead, reducing facial marks while improving stability without hindering your line of sight. Our Advanced Cushion Technology™ (A.C.T.) with AIR^ogel creates a luxurious, ergonomic cushion that reduces pressure on your face.

Move Freely with Our New Full Face Mask



AIR^ogel
Maximum comfort, minimum weight. Our patented AIR^ogel minimizes pressure points and skin irritation.

Advanced Cushion Technology (A.C.T.)
Our ergonomic cushion technology delivers an effective seal for the most comfortable fit.

Quick-Release Headgear
Lightweight, breathable headgear slides off easily with one hand.

360° Rotating Elbow
For optimal tube positioning and easier mobility

Flexible Spacebar™
Unique bendable design provides superior stability, holding a seal without contacting the skin

Size Indicator
Available in three sizes
S M L

Oxygen Port
Allows easy delivery of supplemental oxygen.

Active Headgear Connector
Moves as the patient moves for a reliable seal.



What makes our Innova Full Face Mask different?

- *Flexible Spacebar* is designed to provide superior stability without contacting the skin.
- Our unique *Active Headgear Connector* moves when you move to help maintain a seal throughout the night.
- Available in three sizes with simple and breathable quick-release headgear.

Advanced Cushion Technology (A.C.T.)

Our patented AIR^ogel with an ergonomic design and luxurious feel reduces pressure while maintaining an effective seal.



Active Headgear Connector



Fifteen degrees of rotation provide relief from headgear tension, allowing the mask to move with the patient for a reliable seal.

Innova Full Face Mask with EZfit Headgear Specifications:

Product Number	Case Quantity	Style	Leak Port	Entrainment Valve	Size	Mask Weight	Mask Dimensions
50852	10	Vented	✓	✓	Small	84.3g	3.85"D x 3.90"W x 6.62"H
50853	10	Vented	✓	✓	Medium	89.2g	3.95"D x 4.21"W x 6.75"H
50854	10	Vented	✓	✓	Large	93.0g	4.05"D x 4.36"W x 6.75"H
55051	10	EZFit Headgear	N/A	N/A	One Size	N/A	N/A



Does not contain natural rubber latex

MiniMe²

Sleepnet also offers pediatric masks for our smallest patients



Contact us today.

To place your order, call an authorized distributor.

Call 1.800.742.3646

Visit www.sleepnetmasks.com to learn more.

Sleepnet Corporation, 5 Merrill Industrial Drive,
Hampton, New Hampshire 03842 USA



Scan code to view Innova FF Fitting Tips video.