



CPAP & NIV

Innova[®] Nasal Mask

Innovative fit. Maximum comfort. Freedom of movement

The Innova[®] Nasal Mask, designed with our active headgear connector, is a lightweight nasal mask that provides freedom of movement - improving stability and preventing leaks.



AIR°gel[®]

Sleepnet's AIR°gel[®] is an advanced cushion technology that conforms to the contours of each individual face creating a secure and comfortable seal.

Quick-Release

Our quick release masks feature a slot on the side that allows the user to easily put on or take off the mask and headgear.

360° Rotating Elbow

Designed to provide flexibility and comfort, the 360° rotating elbow allows the user to change sleep positions, move their head, or adjust the mask without discomfort.

Touchless Spacebar


Adding comfort with its contactless design, our headgear attaches to the touchless spacebar and provides additional stability by keeping the mask seal in place when the user moves.

Available Sizes:

S/M M/L

Active Headgear Connector

Ideal for active sleepers, our active headgear connector is designed to move and adjust while maintaining a tight yet comfortable seal and preventing air leaks.





What makes our Innova® Mask different?

Comfort - One of the most significant advantages of AIR°gel® technology is the unparalleled comfort it provides. The gel molds to the user's face, reducing pressure points and maximizing comfort.

Lightweight - Our masks are made from lightweight materials, designed to provide support and structure while ensuring durability and comfort.

Ease of Movement - Our active headgear connector is designed to move and adjust as the wearer moves, allowing the mask to maintain a seal and prevent air leaks. This can be especially helpful for active sleepers, allowing for ease of movement during sleep.

Less Skin Irritation - Our AIR°gel® technology is soft and flexible and helps to reduce irritation.

AIR°gel® Technology:

AIR°gel® is a leading edge innovation developed and designed to enhance the comfort and effectiveness of our respiratory masks. This advanced cushion technology features medical grade silicone gel that conforms to the contours of each individual face creating a secure and comfortable seal that ensures a snug fit while minimizing air leakage around the edges. Our AIR°gel® technology is soft and flexible and helps to reduce irritation, making it suitable for extended wear.



Innova® Nasal Mask with Headgear Specifications:

Product Number	Description	Leak Port	Cushion Size	Mask Weight
50162	Innova® Nasal Vented Mask	✓	S/M	57.0 g
50163	Innova® Nasal Vented Mask	✓	M/L	65.1 g
50592	Innova® Nasal Non-Vented Mask	X	S/M	57.0 g
50593	Innova® Nasal Non-Vented Mask	X	M/L	65.1 g
55051	EZFit Headgear only	N/A	One Size	N/A



Also by Sleepnet®
Mojo® Full Face Mask

Contact us today:

- ☎ 1-800-742-3646
- ✉ info@sleepnetcorp.com
- 🌐 www.sleepnetmasks.com



The following product is protected by US patents 10,974,007, 9,199,054, and other patents pending



FDA
Cleared

Does not contain natural rubber latex

M90008 Version 1

