






PEDIATRIC

MiniMe®2 Nasal Mask

For the little faces that matter most

The MiniMe®2 Nasal Pediatric Mask, developed with guidance from leading children's hospitals, provides a comfortable, non-invasive design that is easy to fit and adjust. It provides a soft, secure, and personalized fit, allowing for unmatched comfort for the extra small faces we care so much about.





What makes our MiniMe®2 different

Softer Comfort - Sleepnet's AIR°gel® is a leading edge innovation developed and designed to enhance the comfort and effectiveness of our respiratory masks.

Personalized fit - Custom Fit Technology allows the user to easily shape the mask by gently bending and molding it along the sides, top, and bottom to provide a consistent and secure fit.

Easy to fit - The mask is specifically designed to provide an optimal fit and superior seal with less pressure - an essential feature in the pediatrics.



Touchless Spacebar:

Our touchless spacebar is a mask component that sits vertically at the top of the mask interface. The headgear attaches to the spacebar without contacting the skin. It is designed to provide additional stability by holding the mask seal in place when a patient moves. It also provides added comfort with its contactless design.

MiniMe®2 Pediatric Mask with Headgear Specifications:

Product Number	Description	Headgear Size	Mask Size	Dead Space	Mask Weight
60216	MiniMe® 2 Pediatric Vented Nasal Mask with Headgear	XS/S	XS	20 ml	51.5 g
60217	MiniMe® 2 Pediatric Vented Nasal Mask with Headgear	XS/S	S	25 ml	49.0 g
60218	MiniMe® 2 Pediatric Vented Nasal Mask with Headgear	M/L	М	37 ml	61.5 g
60221	MiniMe® 2 Pediatric Vented Nasal Mask with Headgear	M/L	L	40 ml	58.5 g
60256	MiniMe® 2 Pediatric Non-Vented Nasal Mask w/Headgear	XS/S	XS	55 ml	48.5 g
60257	MiniMe® 2 Pediatric Non-Vented Nasal Mask w/Headgear	XS/S	S	60 ml	46.0 g
60258	MiniMe® 2 Pediatric Non-Vented Nasal Mask w/Headgear	M/L	М	72 ml	58.5 g
60261	MiniMe® 2 Pediatric Non-Vented Nasal Mask w/Headgear	M/L	L	75 ml	55.0 g
55055	MiniMe® 2 Headgear only	XS/S	N/A	N/A	7.0 g
55056	MiniMe® 2 Headgear only	M/L	N/A	N/A	7.0 g

MiniMe®2 Headgear Sizing Chart:					Recommended Fit				Fits in some cases				
Head Size - Metric	35.5	38	41	43	46	48	51	53.5	56	58.5	61	63.5	
Head Size - Inches	14	15	16	17	18	19	20	21	22	23	24	25	
55055													
55056													



Also by Sleepnet® MiniMe® Pediatric **Contact us today:**



1-800-742-3646



info@sleepnetcorp.com



www.sleepnetmasks.com

Does not contain

natural rubber latex









