



CPAP

Phantom[®]2 Nasal Mask

Your mask, your fit

The Phantom[®]2 Nasal Mask is a moldable, low-profile nasal mask designed for a personalized fit with better visibility for the user. The mask is lightweight, comfortable, and will fit most facial features. The Phantom[®]2 and iQ[®]2 cushions are interchangeable using the same mask frame and headgear. Readily available for resupply.

45° Venting System

Sleepnet's 45-degree venting system is designed to direct air away from the face, and not toward the user or onto their skin during exhalation.

Custom Fit Technology

Sleepnet's Custom Fit Technology is our proprietary technology made from a flexible, pliable material that responds to gentle manipulation and accommodates changes in weight and facial asymmetry.

AIR[°]gel[®]

Sleepnet's AIR[°]gel[®] is an advanced cushion technology that conforms to the contours of each individual face creating a secure and comfortable seal.

360° Rotating Swivel

Designed to provide greater flexibility and enhance comfort, the 360-degree rotating swivel adjusts to the wearer's natural head movement, allowing them to move freely in any direction during sleep.

One Size Cushion

Larger fitting range.


Easy-Release Magnetic Clips

Our magnetic clips are designed to simplify the process of attaching and detaching the mask from the headgear.

Replaceable Cushion

Our replaceable cushions can be removed and cleaned, helping to maintain better hygiene and mask performance.

Available for Resupply



What makes our Phantom[®]2 different


Comfort - Sleepnet's AIR[°]gel[®] is a leading edge innovation developed and designed to enhance the comfort and effectiveness of our respiratory masks.

Personalized fit - Custom Fit Technology allows the user to easily shape the mask by gently bending and molding it along the sides, top, and bottom to provide a consistent and secure fit.

Lightweight - Our masks are made from lightweight materials, designed to provide support and structure while ensuring durability and comfort.

Clear Field of Vision - The iQ[®]2 and Phantom[®]2 masks are designed to offer a Clear Field of Vision giving the user the ability to watch television, read, or wear glasses while the mask stays comfortably in place.

The Phantom[®]2 and iQ[®]2 cushions are interchangeable while using the same mask frame and headgear.



Phantom[®]2 with TruFit headgear and magnetic clips

Product Number	Description	Mask Size	Mask Weight
50182	Phantom [®] 2 Nasal Vented Mask	One Size	75.0 g
50224	Phantom [®] 2 Nasal Cushion	One Size	75.0 g
50204	Swivel Frame / Tubing Assembly	One Size	N/A
55042	TruFit Headgear with Magnetic Clips	One Size	N/A

iQ[®]2 mask

Product Number	Description	Mask Size	Mask Weight
50169	iQ [®] 2 Nasal Vented Mask	One Size	75.0 g
50223	iQ [®] 2 Nasal Cushion	One Size	75.0 g



Also by Sleepnet[®]
iQ[®]2

Contact us today:

- 1-800-742-3646
- info@sleepnetcorp.com
- www.sleepnetmasks.com



The following product is protected by US patents 11,998,694, 10,974,007, 7,938,117, and other patents pending

FDA
Cleared

Does not contain natural rubber latex

M90141 Version 1

